



SCION Thermal Desorber for GC and GC-MS

- Automated, robust and reliable solution for thermal desorption analysis



Introducing the partner of your GC/ GC-MS system for environmental and material analysis

SCION Thermal Desorber coupled with SCION GC/
GC-MS is an ideal platform to measure trace-level
volatile and semi-volatile organic compounds (VOCs
and SVOCs) in various samples: air, flavors,
materials, fragrance, etc.

With a full line of GC and GC-MS instrumentation,
sample introduction devices, custom solutions and
complete selection of consumables and accessories,
SCION is true single-source provider of fully
integrated, world-class GC solutions.


Representative application area of Scion TD-GC/GC-MS system

- TO-1, TO-2, TO-17
- Ambient and Urban Air Monitoring
- Industrial Emissions Monitoring
- Landfill gas analysis
- Odor Analysis
- Photochemical Assessment Monitoring Stations (PAMS)



● Application: EPA Ozone Precursor 56 compounds

Flow diagram of SCION TD-GC



Peak Identification

- | | | | | |
|-------------------|------------------------|----------------------------|----------------------------|----------------------------|
| 1. Ethane | 13. n-Pentane | 25. 2,4-Dimethylpentane | 37. 3-Methylheptane | 49. o-Ethyltoluene |
| 2. Ethylene | 14. trans-2-Pentene | 26. Benzene | 38. n-Octane | 50. 1,2,4-Trimethylbenzene |
| 3. Propane | 15. 1-Pentene | 27. Cyclohexane | 39. Ethylbenzene | 51. n-Decane |
| 4. Propylene | 16. cis-2-Pentene | 28. 2-Methylhexane | 40. m/p-Xylene | 52. 1,2,3-Trimethylbenzene |
| 5. Acetylene | 17. 2,2-Dimethylbutane | 29. 2,3-Dimethylpentane | 41. Styrene | 53. m-Diethylbenzene |
| 6. Isobutane | 18. 2,3-Dimethylbutane | 30. 3-Methylhexane | 42. o-Xylene | 54. p-Diethylbenzene |
| 7. n-Butane | 19. 2-Methylpentane | 31. 2,2,4-Trimethylpentane | 43. n-Nonane | 55. n-Undecane |
| 8. trans-2-Butene | 20. 3-Methylpentane | 32. n-Heptane | 44. Isopropylbenzene | 56. n-Dodecane |
| 9. 1-Butene | 21. Isoprene | 33. Methylcyclohexane | 45. n-Propylbenzene | |
| 10. cis-2-Butene | 22. 1-Hexene | 34. 2,3,4-Trimethylpentane | 46. m-Etyltoluene | |
| 11. Cyclopentane | 23. n-Hexane | 35. Toluene | 47. p-Ethyltoluene | |
| 12. Isopentane | 24. Methylcyclopentane | 36. 2-Methylheptane | 48. 1,3,5-Trimethylbenzene | |



Consumables Highlight – SCION Columns

SCION offers the highest quality, most reproducible GC capillary, PLOT and packed columns, with the most demanding quality assurance specifications in the industry. For GC-MS applications and trace-level analysis, SCION also provide the MS range capillary columns with highly inert and lowest bleed that offer low background, high S/N values and minimal peak tailing.

Accessory and Consumables for SCION TD-GC/GCSQ Analysis

Tube conditioner

Tube conditioner is an ideal device to clean and re-condition the sorbent tube for the use of thermal desorption analysis. Baseline is guaranteed after the re-conditioning and ensuring no carry-over.



Available Various types of Sorbent Tubes

- Tenax-TA
- Tenax-GR
- Carbosieve S-III
- Hayesep porous polymers
- Carbopack, Carboxen, etc.

Focusing trap

- Focusing trap with various sorbent materials
- Depending on polarity, other properties or boiling point range of compounds to be determined various focusing traps are available



For research use only. Not for use in diagnostic procedures.



www.ScionInstruments.com