

Automation of Sample Preparation

SPE-03 8-Channel System


Purpose of Sample Preparation

- Remove Interference from instrumental analysis
- Enrich analytes for lower detection limit
- Convert analytes to a form suitable for analysis

Solid Phase Extraction

- Common method for extraction and clean up
- Much lower solvent usage than liquid/liquid extraction
- Cleaner extracts, less matrix effects

Basic steps of Solid Phase Extraction



Challenges of Manual SPE

Time Consuming – requires constant supervision

Tedious – many steps for conditioning/loading/elution

Challenging – maintaining flow rates, SPE cartridge clogging

Low throughput – 1-2 batches/shift, 4-12 samples/batch

Human Error – inconsistent results, mistakes during extraction



Challenges of Manual SPE

Time Consuming – requires constant supervision

Tedious – many steps for conditioning/loading/elution

Challenging – maintaining flow rates, SPE cartridge clogging

Low throughput – 1-2 batches/shift, 4-12 samples/batch

Human Error – inconsistent results, mistakes during extraction



Automated SPE

Push start and walk away

All conditioning/loading/elution/bottle rinsing steps automated

Constant flow rates, much higher resistance to clogging

Up to 4 batches/shift, 8 samples/batch with 1 system

Consistent results and extraction time, no mistakes



SPE-03 vs Other SPE Systems



SPE-03	Other Systems
8-Channel	1 to 6 Channels
Patented multi-channel valve technology – two valves to handle all 8 samples	Conventional valves – three 3-way valves required for each channel, 24 valves for 8 samples
Compact footprint (34cm x 34cm x 45cm)	Much larger footprint, cannot fit fume hood
Up to 2 fractions per sample	Usually 1 fraction
All steps are performed in parallel	Some only perform sample loading in parallel, conditioning/elution in series
Automated sample bottle rinsing	Requires system pause and manual rinsing
Built-in touch screen computer – up to 100 methods, 100 lines per method, exportable	Requires external computer
Compatible with 1/3/6mL SPE cartridges (other sizes customizable)	Requires separate adapters for 1/3/6mL SPE cartridges

Working Principle

flexible • compact • isolated

Patented Valve Design to Achieve Complex Liquid Handling

PromoChrom's multi-functional valve is based on our flow-path-integration technique. The function of one such valve is equivalent to several normal stream selection valves and isolation valves, keeping the design compact and easy to maintain.




Valve 1 - Stream-selection valve

No contact with samples. Select from 6 solvents, nitrogen, air and bypass waste

Valve 2 - Multi-functional valve

8 Channels, each with dedicated ports for samples, syringe pumps and SPE cartridges



User Interface

computer-free • intuitive • simple

- Step 1. Choose method
- Step 2. Press start

Access different methods

Settings and direct control

Start/Stop/Pause


Select which samples to run


View and edit all method steps directly


Method name and description

setup 3

▼








SPE-03P v1.1


X

Samples

☒ 1 to 4

☒ 5 to 8





Action	Inlet 1	Inlet 2 (ratio)	Flow	Volume
Elute W2	Solvent 1	-	30	3.0
Elute W2	Solvent 2	-	30	3.0
Elute W2	Solvent 3	-	30	3.0
Elute W2	Solvent 4	-	30	3.0
Elute W2	Solvent 5	-	30	3.0
Elute W2	Solvent 6	-	30	3.0
(Add New)				

Method : setup 3 - priming the solvent lines

Water

- Drinking, Ground, Surface, Waste

PFAS, Pesticides, Drugs, PAHs, PCBs



Marine Diesel Oil, Jet Fuel

Saturated hydrocarbon, aromatic hydrocarbon, polar components



Food/Feed

Mycotoxins, pesticides



Extract Clean Up

- Solid, Biosolid and Soil samples

PFAS, EPHs, PPCPs



Protein Purification

Antibodies, antigens, metabolites



Applications

Application Example – PFAS in Drinking Water

EPA Method 537.1

EPA Method 533

DOD

ISO 21675

Modified & Proprietary Methods



- MOD-005 Minimal-Teflon configuration achieving low background
- MOD-004 Automated sample bottle rinsing – a step required by most PFAS methods
- Up to 6 solvents – enough to perform any PFAS method
- Used by the largest commercial labs, government labs and universities in USA
- Can also perform clean up after extraction of PFAS from solids, biosolids and soil samples


Application Example – EPA Method 537.1

Method Details

- Looking at 18 PFAS compounds in drinking water
- 250mL samples extracted using 6mL 500mg SDVB cartridges
- Extracts evaporated to dryness and re-constituted to 1mL 96:4 MeOH/H₂O
- Run extracts on LC-MSMS
- System background requirement **<1/3 MRL**
- System recovery requirement **70% - 130%** for IDC and **50% - 150%** for MRL validation

Extraction Steps

- Precondition with 15mL MeOH, followed by 15mL H₂O
- Add 3mL H₂O
- Load sample at 10-15mL/min
- Rinse sample bottles with 2 aliquots of 7.5mL H₂O and deliver rinsate through SPE cartridges
- Dry SPE cartridges for 5mins
- Rinse sample bottles with 2 aliquots of 4mL MeOH and elute SPE cartridges
- Extraction time of 8 samples takes ~75 minutes on the SPE-03 system



Application Example – EPA Method 537.1

System Background

- Most labs in USA have 2ppt MRL
- Background requirement is $<1/3$ MRL, ie. < 0.667 ppt, all compounds fall well below this limit

Initial Demonstration of Low System Background - 8 x LRBs

Requirements: For 2ppt MRL, background must be < 0.667 ppt

Compound	Pos 1	Pos 2	Pos 3	Pos 4	Pos 5	Pos 6	Pos7	Pos 8
PFBS	ND	ND	ND	ND	ND	ND	ND	ND
PFHxA	ND	ND	ND	ND	ND	ND	ND	ND
HFPO-DA (GenX)	ND	ND	ND	ND	ND	ND	ND	ND
PFHpA	ND	ND	ND	ND	ND	ND	ND	ND
PFHxS	ND	ND	ND	ND	ND	ND	ND	ND
ADONA	ND	ND	ND	ND	ND	ND	ND	ND
PFOA	0.060	0.036	0.052	0.007	0.076	0.072	0.012	0.006
PFOS	ND	ND	ND	ND	ND	ND	ND	ND
PFNA	ND	ND	ND	ND	ND	ND	ND	ND
9Cl-PF3ONS	ND	ND	ND	ND	ND	ND	ND	ND
PFDA	ND	ND	ND	ND	ND	ND	ND	ND
NMeFOSAA	ND	ND	ND	ND	ND	ND	ND	ND
PFUnA	ND	ND	ND	ND	ND	ND	ND	ND
NEtFOSAA	ND	ND	ND	ND	ND	ND	ND	ND
11Cl-PF3OUdSND	ND	ND	ND	ND	ND	ND	ND	ND
PFDoA	ND	ND	ND	ND	ND	ND	ND	ND
PFTTrDA	ND	ND	ND	ND	ND	ND	ND	ND
PFTA	ND	ND	ND	ND	ND	ND	ND	ND

* Results obtained from customer lab during on-site system installation

Application Example – EPA Method 537.1

System Recovery

- IDC requires **70%-130%** and RSD <**20%**
- MRL requires **50%-150%**

IDC - 4 x 50ppt LFBs

Requirements: 70% - 130%, RSD < 20%

Compound	%Recovery	%RSD
PFBS	84	5.16
PFHxA	93	7.81
HFPO-DA (GenX)	95	6.59
PFHpA	104	8.71
PFHxS	99	1.81
ADONA	101	4.92
PFOA	104	5.60
PFOS	95	3.98
PFNA	105	4.73
9CI-PF3ONS	96	1.88
PFDA	96	8.48
NMeFOSAA	101	3.93
PFUnA	96	6.78
NEtFOSAA	101	1.26
11CI-PF3OUdS	86	1.84
PFDaA	87	4.83
PFTTrDA	89	7.81
PFTA	85	10.11

MRL - 7 x 2ppt LFBs

Requirements: 50% - 150%

Compound	%Recovery	%RSD
PFBS	100	5.69
PFHxA	101	4.77
HFPO-DA (GenX)	97	5.14
PFHpA	111	5.02
PFHxS	104	3.79
ADONA	101	5.92
PFOA	112	8.08
PFOS	102	2.24
PFNA	105	7.59
9CI-PF3ONS	96	2.95
PFDA	96	8.52
NMeFOSAA	98	5.56
PFUnA	100	5.47
NEtFOSAA	103	4.21
11CI-PF3OUdS	95	6.38
PFDaA	95	12.66
PFTTrDA	95	11.45
PFTA	92	5.99

* Results obtained from customer lab during on-site system installation

Full SPE-03 Specifications

No. of Samples	8 in parallel
No. of fractions	2
No. of waste channels	2
No. of solvents	6
Sample volume	0.5 – 4000 mL
Fraction volume	Up to 50 mL
SPE cartridge size	1/3/6 mL (Customizable)
Flow rate	0.5 – 100 mL/min
Fluid delivery	Positive pressure
Display	5" resistive touch
No. of methods	100
Method actions	Cartridge pre-condition/soak/wash, add sample, elution, sample bottle rinsing, sample bottle shaking, sample line cleaning, air purge, solvent mixing, nitrogen dry, pause
Dimensions	34 cm x 34 cm x 45 cm
Weight	13 kg
Power	1.5 A @ 24 VDC

Configuration Options

MOD-003

SPE Disk Kit

Expansion for using 47mm disks, option only available with MOD-00P. Includes disk rack and 8 x disk holders.



MOD-003

MOD-004

Sample Bottle Rack

For automated sample bottle rinsing, recommended for up to 250mL bottles. Includes 2 racks and 8 x bottle rinsing adapters.



MOD-004

MOD-005

Minimal-Teflon Option

For PFAS applications. Replaces all PTFE tubing.

MOD-00P

Volume-Matrix Plus Option

Automated sample bottle rinsing with separate sample and rinse lines to handle tough matrices. Recommended for 500mL to 4L samples that require bottle rinsing or sample filtration.



MOD-00P