



Answers for Science.
Knowledge for Life.™




CESI-MS a technique for the characterization and analysis of both intact and digested biopharmaceuticals

Dr. Stephen Lock


Agenda

- The use of CE in Biopharma
- What is CESI-MS?
- Applications using CESI-MS in BioPharma
 - 100% sequence coverage and PTM detection for monoclonal antibody studies.
 - The use of CESI-MS in intact protein analysis in relation to size and charge heterogeneity profiling.
 - How CESI-MS can be used in glycoprotein analysis in the characterisation of EPO & IFN
- Summary


Traditional CE System Schematic



Assessment of IgG Purity and Heterogeneity




1. CE-SDS = Purity analysis



2. CIEF & CZE = charge heterogeneity analysis

“Capillary electrophoresis uses one or more capillaries as migration channels for electrophoresis and increasingly has become the procedure of choice when an electrophoretic separation method is needed. This is because CE is easier to perform, requires less time, and allows better precision and robustness than PAGE.”

USP Guideline for Submitting Requests for Revision to *USP-NF, v3.1*




3. Glycan analysis = microheterogeneity determination

Multi-site Studies in Industry Illustrate Portability of CE

- 2006: CE-SDS Gel: Chromatographia, 64, 359-368; A series of collaborations between various pharmaceutical companies and regulatory authorities concerning the analysis of biomolecules using CE
- 2011: CIEF: Chromatographia, 73, 1137-1144; Intercompany study to evaluate the robustness of C-IEF technology for the analysis of monoclonal antibodies:
- 2012: Imaged CIEF: J. Separation Science, 335, 3124-3129; Robustness of Imaged Capillary IEF methodology for the analysis of monoclonal antibodies: An Inter-laboratory Study
- 2013: N-Glycan Mapping Study: CE Pharm 2013: Pending Publication
- 2015: Evaluation of capillary zone electrophoresis for charge heterogeneity testing of monoclonal antibodies. Journal of Chromatography B, 983–984 (2015) 101–110

Synopsis:


“CE Methods Can Be Reliable, Robust ,and Transportable Across Sites”



What is CESI – MS?

Combining CE & ESI-Mass Spectrometry

CESI 8000 High Performance Separation-ESI Module



Ultra-low flow rates of
< 30nL/min

High resolution,
ultra-low flow rate
CE separations
coupled with
high resolution,
high sensitivity MS




TripleTOF® 5600+
and 6600 Systems


CESI - "The Integration of Capillary Electrophoresis (CE) with Electrospray Ionization (ESI) Into a Single Dynamic Process Within the Same Device"

CE with an Integrated Electrospray Ionization

“ The Integration of Capillary Electrophoresis (CE) With Electrospray Ionization (ESI) Into a Single Dynamic Process Within The Same Device”



Influence of Flow Rate on Sensitivity





Infused sample, Angiotensin I at 2 ng/mL in 10% acetic acid. ESI voltage, -1,350 V; detection range, 50-3000 m/z


Evolution of the peak intensity of Angiotensin I (a) and detection sensitivity as a function of the flow rate (b) Experimental conditions: capillary electrophoresis; bare fused silica capillary with a porous tip, total length 88.5 cm × 30 μm i.d. × 150 μm o.d.; Infused sample, Angiotensin I at 2 ng/mL in 10% acetic acid. Mass spectrometry; capillary voltage, -1350 V; detection range, 50-3000 m/z,

Reducing Ion Suppression Bias at Low Flow Rates

Monitoring maltotetraose suppression in the presence of neurotensin



Analyte suppression decreases logarithmically below 50 nL/min




Multilevel characterization of mAbs by CESI-MS

Multilevel characterization of mAbs by CESI-MS – Bottom up analysis


Evaluation of CESI-MS for the comprehensive characterization of mAb forms

Immunoglobulin G (IgG)




U.S. National Library of Medicine

Peptide mapping
Reduce, alkylate & digest



mAb Characterisation Workflow

- Primary structure characterisation workflow based on bottom-up proteomics strategy



CESI8000 coupled to 5600 TripleTOF MS

Biosimilarity Case Study: Cetuximab

- A single analysis of each sample sufficient to conclude on the complete similarity regarding AA sequence
- Complete sequence coverage is obtained through peptides without miscleavages or PTMs
- CESI-MS/MS enabled to confirm an error, recently reported in the literature

cetuximab

QVQLKQSGPGLVQPSSLSITCTVSGFSLTNYG
VHVVWRQSPGKGLEWLVGIWSSGGNTDYNTPFT
SRLSINKDNSKSKVFFKMNLSQSNDAIYYCAR
ALTYDYEFAYWGQGLTIVTSAASTKGPSVFP
LAPSSKSTSGGTAALGCLVKDYFPEPVTVSWN
SGALTSGVHTFPAVLQSSGLYSLSSVTVPSSS
LGTQTYICNVNHNKPSNTKVDKRVKPKCDKTH
TCPPCPAPPELLGGPSVFLFPPKPKDTLMIS RTP
EVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKT
KPREEQYNSTYRVVSVLTVLHQDWLNGKEYKC
KVSNKALPAPIEKTISKAKGQPREPQVYTLPPSR
EEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPE
NNYKTPPVLDSDGSFFLYSKLTVDKSRWQQGN
VFSCSVMHEALHNHYTQKLSLSLSPGK

DILLTQSPVILSVSPGERVVSFSCRASQSIGTNIH
WYQQRNGSPRLLIKYASESISGIPSRFSGSGS
GTDFTLINSVESEDIADYQCQNNNWPPTFGA
GTKLELKRVAAPSVFIFPPSDEQLKSGTASVV
CLLNNFYPREAKVQWKVDNALQSGNSQESVT
EQDSKDYSLSSSTLTLKADYEKHKVYACEV
THQGLSSPVTKSFNRGEC


 RGAC


cetuximab-B

QVQLKQSGPGLVQPSSLSITCTVSGFSLTNYG
VHVVWRQSPGKGLEWLVGIWSSGGNTDYNTPFT
SRLSINKDNSKSKVFFKMNLSQSNDAIYYCAR
ALTYDYEFAYWGQGLTIVTSAASTKGPSVFP
LAPSSKSTSGGTAALGCLVKDYFPEPVTVSWN
SGALTSGVHTFPAVLQSSGLYSLSSVTVPSSS
LGTQTYICNVNHNKPSNTKVDKRVKPKCDKTH
TCPPCPAPPELLGGPSVFLFPPKPKDTLMIS RTP
EVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKT
KPREEQYNSTYRVVSVLTVLHQDWLNGKEYKC
KVSNKALPAPIEKTISKAKGQPREPQVYTLPPSR
EEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPE
NNYKTPPVLDSDGSFFLYSKLTVDKSRWQQGN
VFSCSVMHEALHNHYTQKLSLSLSPGK

DILLTQSPVILSVSPGERVVSFSCRASQSIGTNIH
WYQQRNGSPRLLIKYASESISGIPSRFSGSGS
GTDFTLINSVESEDIADYQCQNNNWPPTFGA
GTKLELKRVAAPSVFIFPPSDEQLKSGTASVV
CLLNNFYPREAKVQWKVDNALQSGNSQESVT
EQDSKDYSLSSSTLTLKADYEKHKVYACEV
THQGLSSPVTKSFNRGEC


D. Ayoub et al., *mAbs* **2013**, 5, 699-710

Amino Acid Sequence Characterisation



APK
(m/z 315.2039 ; 2+)

CESI-MS detects both very short and long peptides in the same separation



DYFPEPVTVSWNSGALTSGVHTFPAVLQS
SGLYSLSSVVTVPSSSLGTQTYICNVNHKP
SNTK

(63 amino acids ; m/z 1119.898 ; 6+)

Broad Range of Analysis for MAb Glycan ID

Glycopeptides identified as R.EEQYN(Glycan)STYR.V	mAb glycan abbreviation	Glycan mass (Da)	Monoisotopic mass [M+H] ⁺ (Da)	Mass accuracy (ppm)	Average migration time (min) ¹	Average relative abundance (%) ^{1,2}
	FA2G1 (G1F)	1606.5867	2796.0987	-0.41	34.92 ± 0.43	54.49 ± 1.45
	FA2 (G0F)	1444.5339	2634.0459	0.25	34.72 ± 0.42	17.28 ± 0.97
	FA2G2 (G2F)	1768.6395	2958.1515	-0.01	35.19 ± 0.43	8.59 ± 0.40
	A2 (G0)	1298.476	2487.9880	0.67	35.00 ± 0.41	5.58 ± 0.31
	A2G1 (G1)	1460.5288	2650.0408	0.93	34.88 ± 0.43	4.07 ± 0.09
	FA1G1	1403.5073	2591.0037*	3.04	34.45 ± 0.41	1.52 ± 0.06
	FA2G2S1	2059.7349	3249.2469	0.78	34.45 ± 0.42	1.41 ± 0.11
	FA1	1241.4545	2430.9665	0.40	34.76 ± 0.44	1.25 ± 0.06
	M5 (Man5)	1216.4228	2405.9349	0.09	38.18 ± 0.53	1.13 ± 0.66
	FA2G1S1	1897.6821	3085.1785*	3.10	37.93 ± 0.52	1.02 ± 0.48
	A1G1	1257.4494	2446.9614	0.44	38.03 ± 0.48	0.78 ± 0.40
	FA1G1S1	1694.6027	2884.1147	-3.15	34.39 ± 0.42	0.57 ± 0.03
	A2G2 (G2)	1622.5816	2812.0936	-2.77	34.71 ± 0.39	0.55 ± 0.03
	A1	1095.3966	2284.9086	0.57	41.57 ± 0.64	0.51 ± 0.32
	FA2G2S2	2350.8303	3540.3423	1.77	35.11 ± 0.44	0.44 ± 0.05
	A1G1S1	1548.5448	2738.0568	-2.81	37.90 ± 0.51	0.32 ± 0.15
	FA2BG1	1809.6661	2999.1781	-3.00	35.09 ± 0.44	0.17 ± 0.01
	M4A1G1S1	1710.5976	2900.1096	-1.62	34.83 ± 0.43	0.14 ± 0.01
	M5A1	1419.5022	2609.0142	-1.92	37.95 ± 0.52	0.14 ± 0.07
	A2G2S1	1913.677	3103.1890	-1.74	38.05 ± 0.57	0.06 ± 0.03




¹Error was calculated as standard deviations from triplicate technical CESI-MS replicate runs.

²Relative abundances were calculated using peak areas of glycosylated peptides of the same sequence and charge state (-3).

*[M+H]⁺ mass has a -18.0106 Da loss due to water neutral loss.

Glycoform Characterisation: Fc Region

- Fd glycosylation site characterisation



Glycoforms exhibited by the candidate biosimilar are significantly different from cetuximab


➤ 30 % of glycans contains *N*-acetylneuraminic acid




Rejected as biosimilar

PTM Hot Spot Characterization

N-terminal glutamic acid cyclization characterization



Extracted ion electropherograms of peptides HT01 and modified HT01



- CE mechanism separates of peptide with N-terminal glutamic acid cyclization from the unmodified peptide

Results suggest partial
modification of sample



Favorable conditions to estimate
sample modification level


Separation of Aspartic Acid Isomers: Confirmation in Synthetic Peptides

- Aspartic acid isomers separation by CZE confirmed using a synthetic peptide

Synthetic peptides sequences :

NH₂-GLEWIGYISY **D** GTNNYKPSLK-OH


NH₂-GLEWIGYISY **isoD** GTNNYKPSLK-OH



Multilevel characterization of mAbs by CESI-MS – middle down analysis

Immunoglobulin G (IgG)

IgG1, IgG2, & IgG4




U.S. National Library of Medicine


Disulfide mapping
Alkylate & digest
~25 ng of digest analyzed

*Disulfide scrambling between mAb heavy and light chains can result in functional differences.


Level 2



Disulfide linkages identified on IgG4



HC134-LC194




LC88-LC23

LC134-LC194

HC226-HC229


HC261-HC321



Multilevel characterization of mAbs by CESI-MS – top down analysis

Immunoglobulin G (IgG)

IgG1, 2, & 4




U.S. National Library of Medicine



Intact IgG
Charge Heterogeneity
& Reduced Analysis



Level 3



Charge heterogeneity analysis of IgG1 by CESI-MS

CESI-TripleTOF® 6600 MS

- XIC from 150818_N82_IgG1_Run14.wiff...1 to 1360.21 Da, Gaussian smoothed
- XIC from 150818_N82_IgG1_Run14.wiff (...12 to 1431.83 Da, Gaussian smoothed
- XIC from 150818_N82_IgG1_Run14.wiff (...86 to 2339.68 Da), Gaussian smoothed
- XIC from 150818_N82_IgG1_Run14.wiff (...9.5 to 2851.2 Da), Gaussian smoothed
- XIC from 150818_N82_IgG1_Run14.wiff (...91 to 1566.61 Da), Gaussian smoothed
- XIC from 150818_N82_IgG1_Run14.wiff (...7.9 to 3500.0 Da), Gaussian smoothed



pI range – 0.31

***Reverse migration order**


~ 10 mg/mL sample desalted into
50 mM ammonium acetate, pH 4

BGE – 3% acetic acid (tITP-CZE mode)

~3.5 nL injection → ~ 35 ng injected




Charge heterogeneity analysis of IgG1 by CESI-MS

Single analysis at the intact level unifies multiple inferred analyses




Reduced IgG1 CESI-MS analysis

IgG1 light chains



IgG1 heavy chain

EPO Glycoform Heatmap by CESI-MS of Intact Protein



74 distinct glycoforms detected

Sensitive Glycoform Profiling of β -Interferon-1 α (Avonex) and recombinant human Erythropoietin by CESI-TOF-MS – Haselberg et al.



The CESI 8000 Plus – a new range of detection capabilities . . .

LIF



UV/Vis



PDA



CE Assays & Method Development

MS Analysis

Summary

- CE is a robust and reliable technique – well integrated into Biopharma development and control processes
- CESI-MS – Effectively integrates capillary electrophoresis with electrospray ionization serving to reduce/eliminate ion suppression and increase ionization efficiency for:
 - Insulin Impurities
 - Peptide quantitation
 - Multi-level mAb characterization
 - 1) Intact mAb/protein analysis
 - 2) Reduced mAb analysis
 - 3) Peptide mapping and PTM characterization
 - Proteoform/glycoform characterization
- We are just starting to see where this technique will take analytical chemistry.
- Take Home Message - “Flow Matters”

Acknowledgements



Leiden University Medical Center

Andre Deelder
Oleg Mayboroda
Bart Schoenmaker
Anton Heemskeerck
Paul Hensbergen
Oleg Klychnikov
Rico Derks



Innsbruck Medical University

Herbert Lindner
Klaus Faserl
Bettina Sarg
Leopold Kremser
Astrid Devich
Fatma Dikmen
Bernhard Halfinger



Pierre Fabre

Centre d'Immunologie Pierre Fabre

Alain Beck
Elsa Wagner-Rousset
Marie-Claire Janin-Bussat
Daniel Ayoub
Olivier Colas



University of Strasbourg

E. Leize-Wagner
Y.N. François
R. Gahoual
J. Chicher
L. Kuhn
P. Hammann



Utrecht University

Ad de Jong



Beckman Coulter

Jean-Marc Busnel



SCIEX

Bryan Fonslow
Marcia Santos
John Hudson



Trademarks/Licensing

AB Sciex is doing business as SCIEX.

© 2016 AB Sciex. For Research Use Only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.



Any Questions?
stephen.lock@sciex.com

Answers for Science.
Knowledge for Life.™





Answers for Science.
Knowledge for Life.™

