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High-Speed Ion Chromatography

Features of IC-8100 Series


High-throughput analysis within 5 minutes

In combination with a high-speed separation column, measurement of anions and cations can be completed within 5 minutes.


The analysis time can be greatly shortened.

- The measurement time of anions can be shortened from 16 minutes to 5 minutes (compared to our conventional product).
- The measurement time of cations can be shortened from 20 minutes to 5 minutes (compared to our conventional product).

Example of Standard Anions Measurement

Example of Standard Cations Measurement


Eluent and Flow Gradient function

Eluent step gradient and flow step gradient are possible. It is convenient for simultaneous measurement of wide range retention capacity ions. And the cleaning of column can be efficiently performed.

- Analysis time is greatly shortened by eluent step gradient.
- Column cleaning can be efficiently performed by flow step gradient.
- The degasser has two flow lines.


Analysis under eluent gradient condition


Automatic gel replacing suppressor system

Stable and continuous high sensitivity analysis can be achieved by gel suppressor system which consists of suppressor gel and rotary valve.

- Suppressor gel is automatically replaced according to the injection timing of autosampler.
- Regeneration of suppressor is not needed.
- Contaminants from sample do not accumulate in suppressor.
- There is no decrease in sensitivity due to deterioration of the suppressor.
- Low dead volume reduces sample diffusion.
- Suppressor gel low sensor is equipped.
- Two types of suppressor gel bottles are available. (30 mL or 60 mL)


Flow of Suppressor Gel Replacement

Suppressor Gel Bottle (60 mL)

Features of Gel Suppressor System

The gel suppressor system automatically conducts the exchange of suppressor gel according to the injection timing of autosampler. Since the suppressor gel is changed after each measurement, it is not affected by sample contaminants that adsorbed on the suppressor gel. Therefore, it is possible to perform analysis at a stable baseline every time.

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
2.1 Bơm dung dịch:

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Two racks for 50 samples are equipped as standard.
Continuous measurement can be set in the workstation.
Combined with high-throughput analysis, a large number of samples can be processed efficiently in a short time.

- 


Rank 1	Rank 2	Rank 3	Rank 4	Rank 5	Rank 6	Rank 7	Rank 8	Rank 9	Rank 10	Rank 11	Rank 12	Rank 13	Rank 14	Rank 15	Rank 16	Rank 17	Rank 18	Rank 19	Rank 20	Rank 21	Rank 22	Rank 23	Rank 24	Rank 25	Rank 26	Rank 27	Rank 28	Rank 29	Rank 30	Rank 31	Rank 32	Rank 33	Rank 34	Rank 35	Rank 36	Rank 37	Rank 38	Rank 39	Rank 40	Rank 41	Rank 42	Rank 43	Rank 44	Rank 45	Rank 46	Rank 47	Rank 48	Rank 49	Rank 50	Rank 51	Rank 52	Rank 53	Rank 54	Rank 55	Rank 56	Rank 57	Rank 58	Rank 59	Rank 60	Rank 61	Rank 62	Rank 63	Rank 64	Rank 65	Rank 66	Rank 67	Rank 68	Rank 69	Rank 70	Rank 71	Rank 72	Rank 73	Rank 74	Rank 75	Rank 76	Rank 77	Rank 78	Rank 79	Rank 80	Rank 81	Rank 82	Rank 83	Rank 84	Rank 85	Rank 86	Rank 87	Rank 88	Rank 89	Rank 90	Rank 91	Rank 92	Rank 93	Rank 94	Rank 95	Rank 96	Rank 97	Rank 98	Rank 99	Rank 100
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100

[illegible]

50 samples / rack X 2 racks = 100 samples

The autosampler is equipped with automatic dilution function as standard. Automatic dilution analysis is easily set up in workstation.


- Overlay chromatogram of the dilution of standard anion**
-
- The chromatogram displays several peaks. The highest peak is at approximately 2.2 minutes, reaching an absorbance of about 6 AU. Other significant peaks are at approximately 1.8, 2.5, 3.0, and 3.5 minutes. The peaks are labeled with numbers 1 through 6, indicating the dilution series.

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IC-8100-WS

Streamlining routine analysis and variety of functions IC-8100 series dedicated program pursuing "ease of use"

IC-8100-WS workstation is a program with system control, data analysis, and data management functions. Two systems can be controlled, analyzed, and managed with 1 PC.

Easy operation with 3 applications and various tools

- User manager application / user authority setting and registration management
- Acquisition application / system control, data acquisition
- Analysis application / Data analysis, data management
- Log management tool / log confirmation, audit trail tracking
- Validation tool / System evaluation such as noise, detection limit, trueness, precision etc.

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
ES-8100

Product Number: 0024505


Auto eluent generator dedicated for IC-8100EX

The newly developed eluent generator unit stably prepares the eluent and supplies it to the device. Analysis can be done simply by setting concentrated eluent and pure water. Since freshly generated eluent with the constant composition is always supplied, stable analysis can be achieved without fluctuations of retention time through the measurement.

ES-8100

Example of analysis using ES-8100

Reproducibility of SO_4 peak retention time (RT) and peak area

n=60	RT	P. Area	It is possible to supply eluent with constant composition at all time Column : TSKgel SuperIC-Anion HS Conc. Eluent : TSKgel eluent Conc. IC-A-HS-5 Flow Rate : 1.5 mL/min
Average	3.639	28.40	
CV(%)	0.09	0.17	

Changes in retention time (RT) in a long time analysis

RT of SO_4 Peak	min	Retention time does not change even it is used for a long time, thus stable analysis is possible. Column : TSKgel SuperIC-Anion HS Conc. Eluent : TSKgel eluent Conc. IC-A-HS-5 Flow Rate : 1.5 mL/min
After 1 hour	3.643	
After 5 hours	3.639	
After 10 hours	3.635	

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