# VIET NGUYEN

**Technology Service  
Trading Co., Ltd**

Hotline: +84 (0) 817 663300  
Web: [www.vietnguyenco.vn](http://www.vietnguyenco.vn)  
Email: [info@vietnguyenco.vn](mailto:info@vietnguyenco.vn)

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6/30/25



SCIEX TRIPLE Quad

Viet Nguyen Technology Service Trading Co.,Ltd

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**SCIEX - The leader in mass spectrometry solutions**

**VIET NGUYEN**


Viet Nguyen Technology Service Trading Co.,Ltd

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
**NIJ** Technology Transfer Workshop  
National Institute of Justice  
**ESI: Ion Suppression**

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

**ZenoTOF  
7600 system**

Quadrupole LC/MS Market Share, 2021

Enriched Echo® MS System an ultra-high sample throughput system with Acoustic Ejection Mass Spectrometry (AEMS)

The BioPhase 8800 system

Launched the Intabio ZT system, an industry first, to fully integrate icIEF-UV and mass spectrometry

2023

This system was designed to enable new levels of sensitivity for the quantification of previously undetectable compounds

IDL reserpine: 0.14 fg

THE FUTURE

2024

7500+ system  
A new era in robustness and resilience

Explore the system →

Excellence  
Partners

# SCIEX - The leader in mass spectrometry solutions

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**METHOD NUMBER:** C-003.02

**POSTING DATE:** November 14, 2023

**POSTING EXPIRATION DATE:** No expiration date (multi-laboratory validated method).

**PROGRAM AREA:** Mycotoxins

**METHOD TITLE:** Determination of Mycotoxins in Corn, Peanut Butter, and Wheat Flour Using Stable Isotope Dilution Assay (SIDA) and Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS).

### 2023.10. Method

*The use of trade names in this method constitutes neither endorsement nor recommendation by the U. S. Food and Drug Administration (FDA). Citation of instrumentation is based on the original method validation and supporting method verification and platform extension studies.*

A Shimadzu LC-30 system with a binary pump (LC-30AD), autosampler (SIL-30AC) and column oven (CTO-20AC) was used for method development and validation. A Sciex ExionLC and Shimadzu LC-40 have also been tested in method verification.

Mass spectrometric analysis is performed on a **Sciex QTRAP® 4000, Sciex QTRAP® 6500/6500+, or Sciex Triple Quad 6500/6500+** platform in positive ionization mode. Source and compound dependent parameters are summarized below. MS/MS data are collected in a scheduled multiple reaction monitoring (MRM) mode using Analyst® or Sciex OS for instrument control (Sciex). Quantitative data analysis is performed using MultiQuant™ or Sciex OS Version 1.6 or higher (Sciex).

## FOOD CHEMICAL CONTAMINANTS

### Screening of 154 Veterinary Drug Residues in Foods of Animal Origin Using LC-MS/MS: First Action 2020.04

Aurélien Desmarchelier, Thomas Bessaire, Marie-Claude Savoy, Adrienne Tarres, Claudia Mujahid, Andrea Beck, Pascal Mottier and Thierry Delatour\*

#### (s) LC-MS/MS instruments.—

- (1) For Streams A, B, and D, Agilent 1290 SL LC system coupled to a **Sciex 5500 or 6500+** Triple-stage quadrupole mass spectrometer or equivalent. Warning: for Streams A and B, the LC-MS/MS system must also be capable to switch from positive to negative mode within less than 50 ms.

Journal of AOAC INTERNATIONAL, 104(3), 2021, 650–681

doi: 10.1093/jaoacint/qsaa168

Advance Access Publication Date: 22 December 2020

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**Table 5 : Methods mainly used by CVUA Stuttgart**

| Method                                                      | Special remarks on Substances                                                                                                                                                                                                                                                                            | LC-MS/MS                                | Comments                                                                                                                                                                                          |
|-------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>M3: Glyphosate &amp; Co. Hypercarb</b><br>(see 5.6.3)    | Glyphosate<br>AMPA<br>N-Acetyl-AMPA<br>N-Acetyl-Glyphosate<br>Ethephon<br>HEPA<br>Glufosinate<br>N-Acetyl-Glufosinate<br>MPPA<br>Fosetyl-Al<br>Phosphonic acid (screening)<br>Maleic Hydrazide<br>Perchlorate (screening))<br>Chlorate (screening)<br>Cyanuric acid<br>Bialaphos<br>Desmethyl-Dimethoate | Agilent 1200<br><b>Sciex QTRAP 5500</b> | Evaluation via solvent calibration and ILISs except for Bialaphos and N-Acetyl-AMPA (ILIS not yet available)<br>M 1.5 and M 1.6 are currently being tested for their suitability to replace M 1.3 |
| <b>M1.4: Method 1.4 (M 1.4): "Per-ChloPhos" (see 5.6.4)</b> | Perchlorate (quantitative)<br>Chlorate (quantitative)<br>Phosphonic acid (quantitative)<br>Bromide (Screening, quantitative)<br>Bromate (quantitative)<br>Thiocyanate                                                                                                                                    | Agilent 1200<br>Sciex QTRAP 5500        | Mostly employed directly (option: screening by M 1.3, if positive -> M 1.4)<br>Dilution 5-fold<br>Evaluation via solvent calibration and ILISs                                                    |

Please be aware that this old REACH registration data factsheet is no longer maintained. The new ECHA CHEM database has been released by ECHA, and it now contains all information from the transition of ECHA's published data to ECHA CHEM here.

[Access ECHA CHEM](#)

 **ECHA**  
<https://echa.europa.eu>

[ECHA: Homepage](#)

This is unique source of information on the chemicals manufactured and imported in Europe. It covers their hazardous properties, classification and labelling, ...

Use of this information is subject to copyright laws and may require the permission of the copyright holder.

## Di-"isononyl" phthalate

EC number: 249-079-5 | CAS number: 28553-12-0

|                                                                                     |                                               |
|-------------------------------------------------------------------------------------|-----------------------------------------------|
|    | - Toxicological Summary                       |
| <b>General information</b>                                                          | ▶ Toxicokinetics, metabolism and distribution |
|   | ▶ Acute Toxicity                              |
| <b>Classification &amp; Labelling &amp; PBT assessment</b>                          | ▶ Irritation / corrosion                      |
|  | ▶ Sensitisation                               |
| <b>Manufacture, use &amp; exposure</b>                                              | ▶ Repeated dose toxicity                      |
|  | ▶ Genetic toxicity                            |
| <b>Physical &amp; Chemical</b>                                                      | - Carcinogenicity                             |
|                                                                                     | ▶ Toxicity to reproduction                    |
|                                                                                     | ▶ Specific investigations                     |
|                                                                                     | ▶ Exposure related observations in humans     |

**Epider**

Current

Administ

**Admini**

Endpoint:

Type of in

Details on exposure:

Of the 2,582 pregnant women participating in the SELMA study, a first morning urine sample was obtained from 2,356 women (91%) during weeks 9–11 of the pregnancy. These samples were analyzed for the MBP, MEP, MBzP, MEHP, oh-MEHP, oxo-MEHP, cx-MEPP, sum DEHP, oh-MMeOP, oxo-MMeOP, cx-MMeHP, sum DiNP. Briefly, 0.2 mL of urine were added with 0.1 mL of ammonium acetate (pH 6.5) and 0.01 mL glucuronidase (*Escherichia coli*) and thereafter incubated at 37°C for 30 min. Then 0.05 mL was added of a 50:50 (vol:vol) water and acetonitrile solution of labeled (3H or 13C) internal standards of all analyzed compounds. A C18 column (2.1 mm i.d. × 50 mm; Genesis Lightn; Grace, Deerfield, IL, USA) was used before the injector to reduce interference. The phthalate metabolites in the samples were separated on a C18 column (1.5 µm, 2.0 mm i.d. × 30 mm VisionHT; Grace). The mobile phases were water and acetonitrile with 0.08% formic acid. The samples were analyzed on a Shimadzu UFLC system (Shimadzu Corporation, Kyoto, Japan) coupled to a **QTRAP5500 triple quadrupole linear ion trap mass spectrometer equipped with a Turbolon Spray source (AB Sciex, Foster City, CA, USA)**. The samples were analyzed in triplicate, and the mean of the two closest were reported. All samples were analyzed in a randomized order. For quality control of the analyses, chemical blanks and two different in-house prepared quality control samples were analyzed in all sample batches. The limit of detection (LOD) was defined as the concentration corresponding to a peak area ratio of three times the standard deviation of the chemical blanks. The creatinine concentrations were analyzed.

Statistical methods:

Urinary levels of phthalate metabolites were log-transformed to normalize distributions and geometric means (GMs) with 95% confidence interval (CIs).

# SCIEX - The leader in mass spectrometry solutions

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**METHOD NUMBER:** C-003.02

**POSTING DATE:** November 14, 2023

**POSTING EXPIRATION DATE:** No expiration date (multi-laboratory validated method).

**PROGRAM AREA:** Mycotoxins

**METHOD TITLE:** Determination of Mycotoxins in Corn, Peanut Butter, and Wheat Flour Using Stable Isotope Dilution Assay (SIDA) and Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS).

### 2023.10. Method

*The use of trade names in this method constitutes neither endorsement nor recommendation by the U. S. Food and Drug Administration (FDA). Citation of instrumentation is based on the original method validation and supporting method verification and platform extension studies.*

A Shimadzu LC-30 system with a binary pump (LC-30AD), autosampler (SIL-30AC) and column oven (CTO-20AC) was used for method development and validation. A Sciex ExionLC and Shimadzu LC-40 have also been tested in method verification.

Mass spectrometric analysis is performed on a **Sciex QTRAP® 4000, Sciex QTRAP® 6500/6500+, or Sciex Triple Quad 6500/6500+** platform in positive ionization mode. Source and compound dependent parameters are summarized below. MS/MS data are collected in a scheduled multiple reaction monitoring (MRM) mode using Analyst® or Sciex OS for instrument control (Sciex). Quantitative data analysis is performed using MultiQuant™ or Sciex OS Version 1.6 or higher (Sciex).

## FOOD CHEMICAL CONTAMINANTS

### Screening of 154 Veterinary Drug Residues in Foods of Animal Origin Using LC-MS/MS: First Action 2020.04

Aurélien Desmarchelier, Thomas Bessaire, Marie-Claude Savoy, Adrienne Tarres, Claudia Mujahid, Andrea Beck, Pascal Mottier and Thierry Delatour\*

#### (s) LC-MS/MS instruments.—

- (1) For Streams A, B, and D, Agilent 1290 SL LC system coupled to a **Sciex 5500 or 6500+** Triple-stage quadrupole mass spectrometer or equivalent. Warning: for Streams A and B, the LC-MS/MS system must also be capable to switch from positive to negative mode within less than 50 ms.

Journal of AOAC INTERNATIONAL, 104(3), 2021, 650–681

doi: 10.1093/jaoacint/qsaa168

Advance Access Publication Date: 22 December 2020

Article

## THE EUROPEAN FOOD RISK ASSESSMENT FELLOWSHIP PROGRAMME

SERIES 1  
2017-2018

EFSA Journal 2018:16(8):e16083

[www.efsa.europa.eu/efsajournal](http://www.efsa.europa.eu/efsajournal)

### 2.2.2. Part 2

In the second part of the work programme, the fellow was trained to use a robotic workstation Hamilton MICROLAB® STAR equipped with a robotic arm connected to a Multilabel Plate Reader PerkinElmer Victor X3, to develop methods that were used for application of real world sample analysis and further method development and investigation stages. In particular, residual activities of each selected enzyme for the respective OPPs were determined. Assay analysis times were optimised to obtain prime reproducibility of data and to better investigate the effects of unique residual activities. To support the bioenzymatic method and **confirm the robustness**, a conventional mass spectroscopy method was developed and verified using **an AB Sciex QTrap 4,500 triple quadrupole MS system** coupled to a Shimadzu Nextera x2 liquid chromatography system.

# SCIEX - The leader in mass spectrometry solutions

## 2022

### Achieving low parts-per-quadrillion detection limits for PFAS analysis in drinking water

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<sup>2</sup>SCIEX, USA; <sup>3</sup>Phenomenex

#### Abstract

This technical note describes the collaboration between SCIEX and Eurofins Environment Testing to address the 2022 US EPA drinking water health advisory levels (HALs) for 4 per- and polyfluoroalkyl substances (PFAS).

<sup>1</sup> Sample preparation and instrumental methods were developed to achieve low parts-perquadrillion (ppq) HALs—which is pg/L equivalent—for perfluorooctanoic acid (PFOA, 4 ppq), perfluorooctane sulfonic acid (PFOS, 10 ppq), perfluorobutane sulfonic acid (PFBS, 10 ppq), perfluorobutane sulfonamide (PFBSA, 10 ppq), hexafluoropropylene oxide dimer acid (HFPO-DA, or

[Request information →](#)

6/30/25

SCIEX TRIPLE Quad

# SCIEX - The leader in mass spectrometry solutions

# 2022

Our Services About Us Scientific Impact Investors Media Centre Contact Us

Eurofins >> Contact Us >> Worldwide Interactive Map >> USA >> Eurofins Environment Testing – Sacramento

Search >>

Type to search...

Go

## Eurofins Environment Testing – Sacramento

Our Services

Resource library Shop Login Register

Request

Products Applications Training Support Service Stories Community About us

Abstract


Introduction

Key benefits of the method for ultra-

### Achieving low parts-per-quadrillion detection limits for PFAS analysis in drinking water

A collaborative exploration of instrument sensitivity and sample preparation for ultra-trace analysis of PFAS in

<sup>1</sup>Eurofins Environment Testing West Sacra  
<sup>2</sup>SCIEX, USA; <sup>3</sup>Phenomenex

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
Front Chem. 2021; 9: 813359.

Published online 2021 Dec 21. doi: [10.3389/fchem.2021.813359](https://doi.org/10.3389/fchem.2021.813359)

PMCID: PMC8724130

PMID: [34993180](https://pubmed.ncbi.nlm.nih.gov/34993180/)

## Towards Higher Sensitivity of Mass Spectrometry: A Perspective From the Mass

### Analyzers

Chang Li,<sup>1</sup> Shiyong Chu  
and Di Tian<sup>1,\*</sup>


### Abstract

Mass spectrometry  
Recent developmen  
measurement of sul  
made consequently




Introducing a delayed DC ramp to quadrupole can improve the transmission efficiency of the quadrupole ([Brubaker, 1968](#)). A short RF-only quadrupole (pre-filter) is placed in front of the main quadrupole with RF and DC supply. The DC voltage is applied to the main quadrupole only after the ions enter the main quadrupole. Thus, the  $a$  and  $q$  parameters of the targeted ions always remain in the stable region. Losses of ions caused by electric fringe fields at the entrance of the main quadrupole are reduced. The sensitivity is significantly improved. This method is used in many modern commercial mass spectrometers. However, it is not useful when using higher **stability regions**. [Huntley and Reilly \(2020\)](#) simulated the acceptance and transmittance of mass filters driven by **sine and rectangular waveforms**. **When a pre-filter was used or the duty cycle was delayed, the sensitivity of both systems is increased by four times.**

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








# Pushing the limits IonDrive Turbo V – 6500 series



**IonDrive  
Turbo V  
Ion Source**

**QTRAP 6500 System**



**QTRAP 5500 System**

**Turbo V Ion  
Source**

## Pushing the limits E Lens technology for greater sensitivity **7500**

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






① Quantitation via MRM analysis

Enhanced Product Ion (EPI) scan of  $m/z$  343.3 in negative ESI mode

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





**4000  
MRMs**

## XIC of positive and negative ionization

**SCIEX 7500 System**

### ALL QUANTIFIER IONS IN MILLI-Q WATER

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**Guided Optimization - MRM Infusion**

MRM Generator

1. Preparation
2. Set Initial Conditions
3. Optimize Declustering Potential
4. Determine the Product Ions
5. Optimize Collision Energy
6. Optimize Collision Cell Exit Potential
7. Report
8. Select a Method
9. Send MRM to Method Editor

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Guided Optimization - MRM Infusion

MRM Generator

- 1. Preparation
- 2. Set Initial Conditions
- 3. Optimize Declustering Potential
- 4. Determine the Product Ions
- 5. Optimize Collision Energy
- 6. Optimize Collision Cell Exit Potential
- 7. Report
- 8. Select a Method
- 9. Send MRM to Method Editor

Select a Method

Select the mode of MRM method generation

☒ Create a new method

☐ Append MRM transitions to an existing method

The list of Existing Methods contains only those methods that are compatible with the options selected in previous steps.

Existing Method

No compatible methods are available

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#### Project Default Settings

##### Quantitative Processing ▸

##### Qualitative Processing

##### Workspace Layout ▸

Set Project wide defaults for quantitative processing method parameters

##### Method Defaults

Signal to Noise Algorithm Relative Noise ▾

##### Integration Defaults

Integration Algorithm AutoPeak ▾

##### Retention Time (RT)

XIC width 0.02 Da

RT Half Window 30 sec

☐ Optimize RT Half Window

Peak selection by Expected RT ▾

##### Advanced Integration

##### Pre-Processing

☐ Smoothing Low ▾

☒ Noise filter ☐ Moving average

##### Integration

☐ Process by group

Interference resolution 50 %

#CustomersSayNiceThings

“SCIEX OS software  
has allowed us to  
lessen the time spent  
on analyzing data by  
streamlining the  
process.”

Melanie Stauffer, HNL Lab Medicine, United States

SCIEX OS software

☒ Apply to current project

☒ Use for new projects in current data root

More Projects...


SCIEX TRIPLE Quad

Apply

Close


Help


6/30/25

*“over 90% of peaks integrated correctly, in the first pass”*

*“The AutoPeak integration algorithm showed a significant improvement for multiple data types compared to the other algorithms, with over 90% of peaks integrated correctly, in the first pass, without parameter optimization; whereas MultiQuant MQ4 algorithm left a lot of work to be done to optimize the integrations.”*

## EC 657/2002 AOAC – EC 808

| Relative intensity<br>(% of base peak) | EI-GC-MS<br>(relative) | CI-GC-MS, GC-MS <sup>n</sup><br>LC-MS, LC-MS <sup>n</sup><br>(relative) |
|----------------------------------------|------------------------|-------------------------------------------------------------------------|
| > 50 %                                 | ± 10 %                 | ± 20 %                                                                  |
| > 20 % to 50 %                         | ± 15 %                 | ± 25 %                                                                  |
| > 10 % to 20 %                         | ± 20 %                 | ± 30 %                                                                  |
| ≤ 10 %                                 | ± 50 %                 | ± 50 %                                                                  |

## Ion Ratio Flagging Rules Set Up in SCIEX OS

### Variable Tolerance

#### Qualitative Rule

Ion Ratio  
Lower Limit

Ion Ratio  
Upper Limit

 Acceptable  
% Difference  
 Marginal  
% Difference  
 Unacceptable  
% Difference  
 (Tolerance levels based on expected ion ratio calculation)

0

0.1

< 50

< 50

>= 50

0.101

0.2

< 30

< 30

>= 30

0.201

0.5

< 25

< 25

>= 25

0.501

1

< 20

< 20

>= 20

< 20

< 20

>= 20

Custom calculation and filtering packages for SCIEX OS software  
based on EU guidelines

Document N° SANTE/11312/2021

Table 3. Identification requirements for different MS techniques.<sup>2</sup>

| MS detector/Characteristics |                                                                 | Acquisition                                                                                                                                   | Requirements for identification |                                                                                                                               |
|-----------------------------|-----------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------|-------------------------------------------------------------------------------------------------------------------------------|
| Resolution                  | Typical systems (examples)                                      |                                                                                                                                               | minimum number of ions          | additionally                                                                                                                  |
| Unit mass resolution        | Single MS<br>quadrupole, ion trap, TOF                          | full scan, limited m/z range, SIM                                                                                                             | 3 ions                          | S/N ≥ 3 <sup>31</sup><br><br>Analyte peaks from both product ions in the extracted ion chromatograms must fully overlap.      |
|                             | MS/MS<br>triple quadrupole, ion trap, Q-trap, Q-TOF, Q-Orbitrap | selected or multiple reaction monitoring (SRM, MRM), mass resolution for precursor-ion isolation equal to or better than unit mass resolution | 2 product ions                  | Ion ratio from sample extracts should be within <b>±30% (relative)</b> of average of calibration standards from same sequence |

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### Maintenance Kit (LC-40DX3)

Maintenance Kit (LC-40DX3); replaces obsolete part number 5085198.

Part Number: 5085216

**\$828.45**

Quantity

1

### Maintenance Kit (SIL-40C X3)

Maintenance Kit (SIL-40C X3)

Part Number: 5085234

**\$1,827.00**

Quantity

1

### PM Kit (4500/5500/6500)

PM Kit compatible with QTRAP6500, TRIPLE5500, QTRAP6500+, 5500+, Triple3500, TRIPLE4500MD, CITRINE QTRAP, CITRINE QQQ, Triple6500+, QTRAP5500, QTRAP4500MD, TRIPLE6500, TOPAZ instruments

Part Number: 5054105

**\$676.20**

### Oil for Single Stage Roughing Pumps (1 Liter)

This product has replaced part number 1034433; MS30/MS40 - 1 liter (32 fluid oz); Compatible with TRIPLE5500, QTRAP6500+, 4000QTRAP, 7500, API5000, CITRINE QTRAP, API4000, QTRAP5500, TRIPLE5500, Triple6500+, CITRINE QQQ, QTRAP6500, 5500+ instruments

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
**SCIEX OS**  
**Software**

A SIMPLE SWITCH  
FOR A BIG DIFFERENCE

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## Brands distributions by Viet Nguyen Co.,Ltd

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## Certificate of Completion

**Nguyen Khac Manh**

**SCIEX Triple Quad 5500+ LCMS/MS Service Training**

10/15/2023

SCIEX confirms that the above-named engineer has been trained to support installation, maintenance, repair & qualification tasks. Eligible to receive continuing training on installation, repair, maintenance, qualification, and has the ability to train customers on operation of the instruments listed in this record. Through the course of employment with SCIEX or by an SCIEX Authorized Dealer, this engineer has access to the latest SCIEX product information, service tools, and spare parts.

This record is being issued on the date noted above, and remains valid provided that the above-named engineer remains employed by SCIEX or by one of SCIEX authorized dealers. Please feel free to contact [CustomerCare@sciex.com](mailto:CustomerCare@sciex.com) to confirm the ongoing validity of this certificate at any time.

*Joe Lepore*  
Global Director  
Commercial And Service Training

## Authorized Service Provider Certificate of Completion

**Nguyen Khac Manh**

Date Issued: 9/18/2023

**Is Trained to Perform Installation Services**

**Mass Spectrometers: 3500/4500/5500/6500 Models**  
**LC Systems: Shimadzu & Exlon Models**

SCIEX confirms that the above-named engineer has been trained on Installation is eligible to receive continuing training on installation and has the ability to train customers on operation of the instruments listed in this record. Through the course of employment with SCIEX as an SCIEX Authorized Service Provider, this engineer has access to the latest SCIEX product information.

This record is being issued on the date noted above, **and remains valid for one year** provided that the above-named engineer remains employed by SCIEX as an Authorized Service Provider.

Please feel free to contact [CustomerCare@sciex.com](mailto:CustomerCare@sciex.com) to confirm the ongoing validity of this certificate at any time.

*Joe Lepore*  
Global Director  
Commercial And Service Training

Satish Chandra Kumar Guttalapalli  
FSE Manager

## Certificate of Completion

**Tran The Duy**

**SCIEX Triple Quad 5500+ LCMS/MS Service Training**

7/16/2023

SCIEX confirms that the above-named engineer has been trained and is eligible to receive continuing training on installation, repair, maintenance, qualification, and has the ability to train customers on operation of the instruments listed in this record. Through the course of employment with SCIEX or by an SCIEX Authorized Dealer, this engineer has access to the latest SCIEX product information, service tools, and spare parts.

This record is being issued on the date noted above, and remains valid provided that the above-named engineer remains employed by SCIEX or by one of SCIEX authorized dealers. Please feel free to contact [CustomerCare@sciex.com](mailto:CustomerCare@sciex.com) to confirm the ongoing validity of this certificate at any time.

*The Lyngere*

Global Director  
Commercial And Service Training

## Certificate of Completion

**Tran The Duy**

**TripleQuad/QTRAP "BlueLine" MS - Classroom Training**

9/21/2023

SCIEX confirms that the above-named engineer has been trained and is eligible to receive continuing training on installation, repair, maintenance, qualification, and has the ability to train customers on operation of the instruments listed in this record. Through the course of employment with SCIEX or by an SCIEX Authorized Dealer, this engineer has access to the latest SCIEX product information, service tools, and spare parts.

This record is being issued on the date noted above, and remains valid provided that the above-named engineer remains employed by SCIEX or by one of SCIEX authorized dealers. Please feel free to contact [CustomerCare@sciex.com](mailto:CustomerCare@sciex.com) to confirm the ongoing validity of this certificate at any time.

*The Lyngere*

Global Director  
Commercial And Service Training

## Certificate of Completion

**Le Anh Tuan**

**SCIEX Triple Quad 5500+ LCMS/MS Service Training**

10/12/2023

SCIEX confirms that the above-named engineer has been trained to support installation, maintenance, repair & qualification tasks. Eligible to receive continuing training on installation, repair, maintenance, qualification, and has the ability to train customers on operation of the instruments listed in this record. Through the course of employment with SCIEX or by an SCIEX Authorized Dealer, this engineer has access to the latest SCIEX product information, service tools, and spare parts.

This record is being issued on the date noted above, and remains valid provided that the above-named engineer remains employed by SCIEX or by one of SCIEX authorized dealers. Please feel free to contact [CustomerCare@sciex.com](mailto:CustomerCare@sciex.com) to confirm the ongoing validity of this certificate at any time.

## Authorized Service Provider Certificate of Completion

**Le Anh Tuan**

Date Issued: 9/18/2023

**Is Trained to Perform Installation Services**

**Mass Spectrometers: 3500/4500/5500/6500 Models**  
**LC Systems: Shimadzu & Exlon Models**

SCIEX confirms that the above-named engineer has been trained on installation is eligible to receive continuing training on installation and has the ability to train customers on operation of the instruments listed in this record. Through the course of employment with SCIEX as an SCIEX Authorized Service Provider, this engineer has access to the latest SCIEX product information.

This record is being issued on the date noted above, and remains valid for one year provided that the above-named engineer remains employed by SCIEX as an Authorized Service Provider.

Please feel free to contact [CustomerCare@sciex.com](mailto:CustomerCare@sciex.com) to confirm the ongoing validity of this certificate at any time.


Global Director


**Commercial And Service Training**

Satish Chandra Kumar Guttalapalli

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


# Field Service Support Software

Provide an overview and supplemental information for the Field Service Support Software used by SCIEX Field Service Engineers (FSEs) and Channel Partners via a Kanguru™ Flash Drive to configure, test, troubleshoot and maintain the Blueline Series Mass Spectrometers.

- Analyst® Service Diagnostics (ASD)
- Tuning Tools (TT)
- Field Service Data Collector (FSDC)
- Agilent T-plus

## QTrap4500 Mass Spectrometer Diagnostic Report

10/13/2023 3:55 PM

User Account: localadmin

Instrument Model: QTrap4500

Ion Source: TurbolonSpray

Serial Number: BW20141504

Target Instrument: QTrap4500

Optimization Set: QTrap4500 Diagnostic

Doc Num Revision: D5197485 A

| Name                | Minimum | Maximum | Readback | Status |
|---------------------|---------|---------|----------|--------|
| Ambient Temperature | 15.00   | 35.00   | 23.87    | Pass   |

Check the Pressure at CAD 0

Pass

The pressure 6.9E-06 Torr was within tolerance levels

### 3. Low Voltage DC Power Supply

| Step                                 | Result    |         | Details |          |        |
|--------------------------------------|-----------|---------|---------|----------|--------|
| Verify Low Voltage DC Power supplies | Pass      |         |         |          |        |
| Name                                 | Set Point | Minimum | Maximum | Readback | Status |
| Electronics DC Power Supply +19.0    | N/A       | 17.00   | 21.00   | 19.19    | Pass   |
| Electronics DC Power Supply +24.0    | N/A       | 22.00   | 26.00   | 23.81    | Pass   |
| Electronics DC Power Supply +3.3     | N/A       | 3.10    | 3.50    | 3.30     | Pass   |
| Electronics DC Power Supply +5.0     | N/A       | 4.50    | 5.50    | 5.00     | Pass   |
| Electronics DC Power                 | N/A       | -21.00  | -17.00  | -18.90   | Pass   |

Disclaimer & copyright:

- This video is intended for MS Vision trained & certified customers only.
- Please observe applicable safety precautions; make sure the instrument is vented, the power is switched off, and disconnected from the mains. Please wear gloves when needed.
- Spectrometry Vision B.V. cannot be held liable for any damages or personal injuries that may occur.
- All images & video content are sole Copyright of Spectrometry Vision B.V. and in such may not be downloaded, reproduced, copied, projected, defaced, edited or used in any way without written permission of Spectrometry Vision B.V.

## Việt Nguyễn – vietCALIB® PREMIUM OF AFTER SALES SERVICE

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
- Increased coverage from limited samples
- High-resolution separation of polar analytes and critical PTMs
- In-depth characterization of proteoforms and protein complexes

**See What You've Been Missing**

PLUS HIGH PERFORMANCE SEPARATION - ESI MODULE

Perfluorocarbons  
PFCs

Pesticide

## On April 10, 2024, EPA announced the final National Primary Drinking Water Regulation (NPDWR) for six PFAS.

EPA established legally enforceable levels, called Maximum Contaminant Levels (MCLs), for six PFAS in drinking water: PFOA, PFOS, PFHxS, PFNA, and HFPO-DA as contaminants with individual MCLs, and PFAS mixtures containing at least two or more of PFHxS, PFNA, HFPO-DA, and PFBS using a Hazard Index MCL to account for the combined and co-occurring levels of these PFAS in drinking water. EPA also finalized health-based, non-enforceable Maximum Contaminant Level Goals (MCLGs) for these PFAS.

| Compound                                                          | Final MCLG                   | Final MCL (enforceable levels)                        |
|-------------------------------------------------------------------|------------------------------|-------------------------------------------------------|
| PFOA                                                              | Zero                         | 4.0 parts per trillion (ppt) (also expressed as ng/L) |
| PFOS                                                              | Zero                         | 4.0 ppt                                               |
| PFHxS                                                             | 10 ppt                       | 10 ppt                                                |
| PFNA                                                              | 10 ppt                       | 10 ppt                                                |
| HFPO-DA (commonly known as GenX Chemicals)                        | 10 ppt                       | 10 ppt                                                |
| Mixtures containing two or more of PFHxS, PFNA, HFPO-DA, and PFBS | 1 (unitless)<br>Hazard Index | 1 (unitless)<br>Hazard Index                          |

For Research Use Only. Not for Use In Diagnostic Procedures.

Food and Environmental

## Analysis of PFAS at low ppt levels in drinking water via EPA method 533

*Using the SCIEX Triple Quad 4500 system*

Matthew Standland, Craig M. Butt, and Matthew Noestheden  
SCIEX, USA

Per- and polyfluorinated alkyl substances (PFAS) are a group of manmade chemicals that have been used for decades in a large host of applications. A recent review focused on 1,400 unique PFAS species<sup>1</sup> although current OECD databases have upwards of 4,500 unique PFAS compounds.<sup>2</sup> The resilience of some PFAS to degradation and their potential for biological harm have spurred health and environmental concerns amongst regulatory agencies. The CDC's national biomonitoring program has found evidence of widespread presence of certain PFAS residues in human serum<sup>3</sup>. One possible route of human exposure is from PFAS contaminated water.<sup>4</sup>

PFAS are known to enter the water supply and bioaccumulate in watershed ecosystems. This can be a health issue for humans who directly consume contaminated water. The US EPA published the fifth unregulated contaminant monitoring rule

(UCMR5)<sup>5</sup>, which proposes maximum residue limits (MRL) for PFAS residues in EPA method 533. Representative chemical structures are shown in Figure 1.


- EPA method 533. 6
- 250 mL of water was passed through Phenomenex Strata-X-AW Weak Anion 500 mg/6 mL.
- Eluted with methanol containing 2% ammonium hydroxide, dried down under nitrogen gas, and reconstituted in 1 mL of 80:20 methanol/ water

**Làm giàu 250 lần**

Chromatography: ExionLC AC system using a 12-minute gradient.


- The Phenomenex Gemini C18 (100 x 3.0 mm, 3.0  $\mu$ m, 00D-4439-Y0)
- The delay column used was a Phenomenex Luna C18 (30 x 3.0 mm, 00A-4252-Y0).
- **The injection volume was 5  $\mu$ L**
- Column temperature of 40°C using a flow rate of 0.6 mL/min.
- Mobile phases consisted of 10 mM ammonium acetate and methanol with 10 mM ammonium acetate.

**Figure 3. TIC chromatogram of 15 ng/mL solvent standard showing chromatographic separation of EPA method 533 PFAS analytes.** Represent target residues are displayed here as extracted from their corresponding MRM in the mass spectrometry method.

**Table 1. Accuracy and precision of method.** Precision CV measurements were based on the area response of the quantifier transition obtained from 7 replicates of water fortified at 6 ng/L. Mean % area values are based on 7 replicates of water fortified at 4 ng/L. Verified MRLs were established using EPA method 533 criteria. All residues passed data quality requirements in section 9 of EPA method 533.

| Compound     | Verified MRL (ng/L) | CV (%) | Mean accuracy (%) | UCMR5 MRL (ng/L) |
|--------------|---------------------|--------|-------------------|------------------|
| 11CI-PF3OUdS | 2                   | 6.6    | 85.2              | 5                |
| 9CI-PF3ONS   | 2                   | 7.4    | 96.8              | 2                |
| ADONA        | 2                   | 3.5    | 94.9              | 3                |
| HFPO-DA      | 2                   | 5.4    | 94.5              | 5                |
| NFDHA        | 2                   | 4.3    | 91.1              | 20               |
| PFBA         | 2                   | 7.2    | 109.8             | 5                |
| PFBS         | 2                   | 6.6    | 97.9              | 3                |
| 8:2 FTS      | 3                   | 9.5    | 117.0             | 5                |
| PFDA         | 2                   | 4.5    | 96.2              | 3                |
| PFDoA        | 2                   | 6.1    | 102.6             | 3                |
| PFEESA       | 2                   | 7.3    | 102.6             | 3                |
| PFHpS        | 2                   | 7.5    | 98.8              | 3                |
| PFHpA        | 2                   | 7.5    | 104.4             | 3                |
| 4:2 FTS      | 3                   | 6.5    | 122.2             | 3                |
| PFHxS        | 2                   | 7.7    | 99.9              | 3                |
| PFHxA        | 2                   | 4.2    | 98.4              | 3                |

**Table 1. Accuracy and precision of method.** Precision CV measurements were based on the area response of the quantifier transition obtained from 7 replicates of water fortified at 6 ng/L. Mean % area values are based on 7 replicates of water fortified at 4 ng/L. Verified MRLs were established using EPA method 533 criteria. All residues passed data quality requirements in section 9 of EPA method 533.

| Compound     | Verified MRL (ng/L) | CV (%) | Mean accuracy (%) | UCMR5 MRL (ng/L) |
|--------------|---------------------|--------|-------------------|------------------|
| 11CI-PF3OUdS | 2                   | 6.6    | 85.2              | 5                |
| 9CI-PF3ONS   | 2                   | 7.4    | 96.8              | 2                |
| ADONA        | 2                   | 3.5    | 94.9              | 3                |
| HFPO-DA      | 2                   | 5.4    | 94.5              | 5                |
| NFDHA        | 2                   | 4.3    | 91.1              | 20               |
| PFBA         | 2                   | 7.2    | 109.8             | 5                |
| PFBS         | 2                   | 6.6    | 97.9              | 3                |
| 8:2 FTS      | 3                   | 9.5    | 117.0             | 5                |
| PFDA         | 2                   | 4.5    | 96.2              | 3                |
| PFDoA        | 2                   | 6.1    | 102.6             | 3                |
| PFEESA       | 2                   | 7.3    | 102.6             | 3                |
| PFHpS        | 2                   | 7.5    | 98.8              | 3                |
| PFHpA        | 2                   | 7.5    | 104.4             | 3                |
| 4:2 FTS      | 3                   | 6.5    | 122.2             | 3                |

| Compound | Verified MRL (ng/L) | CV (%) | Mean accuracy (%) | UCMR5 MRL (ng/L) |
|----------|---------------------|--------|-------------------|------------------|
| PFHxS    | 2                   | 7.7    | 99.9              | 3                |
| PFHxA    | 2                   | 4.2    | 98.4              | 3                |
| PFMPA    | 2                   | 3.4    | 100.4             | 3                |
| PFNA     | 2                   | 4.2    | 93.3              | 4                |
| 6:2 FTS  | 3                   | 5.1    | 112.9             | 5                |
| PFOS     | 2                   | 8.9    | 98.8              | 4                |
| PFOA     | 2                   | 8.0    | 100.0             | 4                |
| PFPeA    | 2                   | 4.2    | 101.4             | 3                |
| PFPeS    | 2                   | 7.1    | 98.5              | 4                |
| PFUnA    | 2                   | 4.1    | 100.2             | 2                |

For Research Use Only. Not for Use In Diagnostic Procedures.

## Food and Environmental

## Quantitation of PFASs in water samples using LC-MS/MS large-volume direct injection and solid phase extraction


Simon Roberts<sup>1</sup>, KC Hyland<sup>1</sup>, Craig Butt<sup>1</sup>, Scott Krepich<sup>2</sup>, Eric Redman<sup>3</sup>, and Christopher Borton<sup>1</sup>

<sup>1</sup>SCIEX, USA; <sup>2</sup>Phenomenex, USA; <sup>3</sup>TestAmerica Laboratories, Sacramento, USA

PFASs are unique chemicals whose physicochemical properties make them important for use in a variety of industrial and consumer products including carpets, cookware, food packaging, fire suppressants, and others<sup>1</sup>. Chemically, PFASs are aliphatic structures containing one or more C atoms on which H substituents have been replaced by F atoms. Classification and naming is typically by the particular functional group present, such as carboxylic acids, sulfonates, phosphonic acids, etc., as well as the length of the carbon chain. Desirable in various industrial applications for their chemical stability and low reactivity, these properties also make PFASs highly resistant to degradation in aquatic environments. Typical concentrations of PFASs found in various environmental water sources range from pg/L to µg/L levels<sup>2</sup>.

This application note presents the methods for the quantitation

**Figure 4. Method 2 chromatography.** Chromatogram of a 10 ng/L matrix spike into groundwater that was diluted with methanol and injected according to Method 2.

# PFAS/PFOS SCIEX 5500 System

**Table 5. Calibration curves for method 1 and 2.** Sensitivity and linearity from 25 to 20,000 ng/L and 1 to 200 ng/L (coefficient of regression, R2) using Method 1 and Method 2, respectively. S/N calculated using MultiQuant Software 3.0.2.

| Compound      | Method 1                 |                         |                         |                              | Method 2                 |                         |                        |                             |
|---------------|--------------------------|-------------------------|-------------------------|------------------------------|--------------------------|-------------------------|------------------------|-----------------------------|
|               | Calibration range (ng/L) | Linear correlation (R2) | S/N of 25 ng/L standard | Accuracy of 25 ng/L standard | Calibration range (ng/L) | Linear correlation (R2) | S/N of 1 ng/L standard | Accuracy of 1 ng/L standard |
| <b>PFCA</b> s |                          |                         |                         |                              |                          |                         |                        |                             |
| PFBA          | 25-20,000                | 0.997                   | 108                     | 104%                         | 1-200                    | 0.997                   | 328                    | 97%                         |
| PFPeA         | 25-20,000                | 0.998                   | 88                      | 103%                         | 1-200                    | 0.999                   | 137                    | 101%                        |
| PFHxA         | 25-20,000                | 0.998                   | 104                     | 93%                          | 1-200                    | 0.999                   | 284                    | 101%                        |
| PFHpA         | 50-20,000                | 0.999                   | 116                     | 101%                         | 1-200                    | 0.993                   | 267                    | 96%                         |
| PFOA          | 25-20,000                | 0.999                   | 117                     | 106%                         | 1-200                    | 0.999                   | 113                    | 99%                         |
| PFNA          | 25-20,000                | 0.990                   | 91                      | 109%                         | 1-200                    | 0.999                   | 137                    | 101%                        |
| PFDA          | 25-20,000                | 0.998                   | 103                     | 105%                         | 1-200                    | 0.997                   | 176                    | 96%                         |
| PFUdA         | 25-20,000                | 0.995                   | 84                      | 101%                         | 1-200                    | 0.998                   | 168                    | 99%                         |
| PFDoA         | 25-20,000                | 0.998                   | 60                      | 101%                         | 1-200                    | 0.994                   | 127                    | 94%                         |
| PFTTrDA       | 25-20,000                | 0.998                   | 32                      | 104%                         | 1-200                    | 0.995                   | 125                    | 95%                         |
| PFTeDA        | 25-20,000                | 0.994                   | 15                      | 107%                         | 1-200                    | 0.998                   | 56                     | 98%                         |
| PFHxDA        | 25-20,000                | 0.999                   | 21                      | 103%                         |                          |                         |                        |                             |
| PFODA         | 25-20,000                | 0.999                   | 33                      | 102%                         |                          |                         |                        |                             |
| <b>PFSA</b> s |                          |                         |                         |                              |                          |                         |                        |                             |
| PFBS          | 25-20,000                | 0.995                   | 31                      | 92%                          | 2-200                    | 0.994                   | 1178                   | 100%                        |
| PFHxS         | 25-20,000                | 0.999                   | 604                     | 103%                         | 1-200                    | 0.998                   | 229                    | 96%                         |
| PFHpS         | 25-20,000                | 0.997                   | 103                     | 105%                         | 1-200                    | 0.999                   | 327                    | 99%                         |
| PFOS          | 25-20,000                | 0.995                   | 312                     | 105%                         | 1-200                    | 0.999                   | 251                    | 99%                         |
| PFDS          | 25-20,000                | 0.998                   | 88                      | 102%                         | 1-200                    | 0.999                   | 516                    | 98%                         |

## An ultra-high sensitivity analysis of PFAS compounds in multiple water sources

### *Using the SCIEX 7500 System*

**Bertram Nield<sup>1</sup>, Jack Steed<sup>2</sup>, Rens-Jan Jaspers<sup>3</sup>, Kay Hup<sup>3</sup>, Abdessamad Chahbouni<sup>3</sup>, Jianru Stahl-Zeng<sup>4</sup>**

<sup>1</sup>SCIEX, The Netherlands; <sup>2</sup>SCIEX, UK; <sup>3</sup>Het Waterlaboratorium; <sup>4</sup>SCIEX, Germany

In this method sub-parts per trillion (ppt) levels of detection for 26 PFAS compounds was achieved with LOD values of 0.2 ng/L in diluent and calculated method detection limits in drinking, ground, and surface water ranging between 0.06 ng/L and 1.12 ng/L.

PFAS compounds are ubiquitous in our environment due to overuse and their lack of breakdown, ensuring that this will be a challenge for decades to come.<sup>1</sup> Therefore, it is imperative to provide rigorous and sensitive analytical testing to regulate these compounds and try to limit their possible effects on human health.

In December 2020, the European Parliament and Council of the European Union released a new directive that sets the limit of PFAS in drinking water to 0.5 µg/L for all PFAS compounds.

### Key features of PFAS analysis using the SCIEX 7500 System

- Sample preparation: **2.5 mL of the water sample + 2 mL of a 50:50 acetonitrile/methanol + 0.225 % formic acid**
- Chromatography: Phenomenex Luna Omega PS C18, 100 Å, 100 x 2.1 mm, 3 µm (PN: 00D-4758-AN),  
Phenomenex Gemini C18, 110 Å, 100 x 2.0 mm, 3 µm (PN: 00D-4439-B0) delay column.
- Mobile phase A was 20 mM ammonium acetate in water and mobile phase B was methanol.
- Mass spectrometry: The analysis was performed using the SCIEX 7500 System, operated in electrospray ionization in negative ion mode.


Figure 1. Extracted ion chromatograms (XICs) of three PFAS compounds at their respective LOQs. From left to right: PFBS, PFHxA and PFOS - all at 0.2 ng/L. The image above shows excellent sensitivity has been achieved with this method at low ng/L levels.


Figure 4. Calibration curve for PFNA. The image above shows the linear range of 0.5 – 50 ng/L for PFNA, highlighting an r value of 0.99990.

Table 1. LOQ of each compound analyzed in diluent with peak-to-peak S/N added to show that each peak is quantifiable, with S/N above 10 x.

| Compound name | LOQ (ng/L) | Peak-to-peak S/N |
|---------------|------------|------------------|
| PFHxS         | 0.20       | 34.5             |
| PFNA          | 0.50       | 25.5             |
| PFOS          | 0.20       | 34.5             |
| PFOA          | 0.50       | 14.3             |
| PFBA          | 0.20       | 11.7             |
| PFPeA         | 0.20       | 19.6             |
| PFHxA         | 0.20       | 24.7             |
| PFHpA         | 0.20       | 23.6             |
| PFDA          | 0.50       | 15.1             |
| PFUdA         | 0.50       | 14.6             |
| PFDoA         | 0.50       | 11.1             |
| PFTTrDA       | 1.00       | 55.4             |
| HFPO-DA       | 0.20       | 31.2             |

## Achieving low parts-per-quadrillion detection limits for PFAS analysis in drinking water

*A collaborative exploration of instrument sensitivity and sample preparation for ultra-trace analysis of PFAS on the SCIEX 7500 system*

Thép Phomsopha<sup>1</sup>, Simon Roberts<sup>2</sup>, Craig M. Butt<sup>2</sup>, Karl Oetjen<sup>2</sup>, Matt Noestheden<sup>2</sup>, Sam Lodge<sup>3</sup>, Andrew N. Patterson<sup>1</sup>

<sup>1</sup>Eurofins Environment Testing West Sacramento;

<sup>2</sup>SCIEX, USA; <sup>3</sup>Phenomenex

This technical note describes the collaboration between SCIEX and Eurofins Environment Testing to address the 2022 US EPA drinking water health advisory levels (HALs) for 4 per- and polyfluoroalkyl substances (PFAS).<sup>1</sup> Sample preparation and instrumental methods were developed to achieve low parts-per-quadrillion (ppq) HALs—which is pg/L equivalent—for perfluorooctanoic acid (PFOA, 4 ppq), perfluorooctane sulfonic acid (PFOS, 20 ppq), perfluorobutane sulfonic acid (PFBS, 2 ppb) and hexafluoropropylene oxide dimer acid (HFPO-DA, or GenX, 10 ppt). The high sensitivity of the SCIEX 7500 system

EPA drinking water HALs demand an unprecedented level of cleanliness and instrumental sensitivity which ultimately necessitates newer, robust analysis techniques.


**Key benefits of the method for ultra-trace level PFAS analysis**

- Sample preparation: The sample extraction followed a modified EPA Method 533 procedure
- 250 mL of sample., load on (SPE) cartridges were the Phenomenex Strata™ X-AW (500 mg).
- The final SPE eluent was adjusted to 10 mL of 80:20 methanol/water (0.1% ammonium hydroxide) and an aliquot was transferred to a polypropylene vial for analysis.

Chromatography: SCIEX ExionLC system

- The analytical column was the Gemini C18 (3  $\mu$ m, 100 x 3 mm, Phenomenex, P/N: 00D-4439-Y0) with a Zorbax Diol guard cartridge (6  $\mu$ m, 12.5 x 4.6 mm, P/N: 820950-911).
- The delay column was the Luna Omega PS C18 (5  $\mu$ m, 50 x 4.6 mm, Phenomenex, P/N: 00B4753-E0).
- The mobile phases were water and methanol (modified with 5mM and 2mM ammonium acetate, respectively) with a flow rate of 0.7 mL/min. The column oven was 45 oC,
- The injection volume was 100  $\mu$ L.

# PFAS/PFOS SCIEX 7500 System


Figure 2. XIC chromatograms for HFPO-DA ( $m/z$  285.0>169.0), PFBS ( $m/z$  298.9>80.0) and PFOS ( $m/z$  499>99.0) quantifier MRM transitions from method extraction samples. Instrument blank, extraction blank, and two spike levels are shown.

**U.S. FOOD & DRUG  
ADMINISTRATION**

## FDA Foods Program Compendium of Analytical Laboratory Methods: Chemical Analytical Manual (CAM)

**METHOD NUMBER:** C-010.03

**POSTING DATE:** 04/12/2024

**POSTING EXPIRATION DATE:** 04/12/2026

**PROGRAM AREA:** Chemical Contaminants

**METHOD TITLE:** Determination of 30 Per and Polyfluoroalkyl Substances (PFAS) in Food and Feed using Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS)

**VALIDATION STATUS:** Single Laboratory Validation (L2) per the [Guidelines for the Validation of Chemical Methods for the FDA Foods Program 3<sup>rd</sup> Edition](#)

### 2021.9 CALCULATIONS

Example calculation for concentration measured on LC-MS/MS to concentration in 5 grams of food with a final extract of 10 mL:

- The lowest calibration curve point is 0.01 ng/mL in 0.5 mL of solution.

$$\frac{0.01 \text{ ng}}{\text{mL}} * 0.5 \text{ mL} = 0.005 \text{ ng}$$

- Since there are 0.005 ng in 0.5 mL of extract, there would be 0.1 ng in the total 10 mL extract

$$0.005 \text{ ng} * \frac{10 \text{ mL}}{0.5 \text{ mL}} = 0.1 \text{ ng}$$

- With a 5 gram food sample, this is equivalent to 20 ng/kg in foods

$$\frac{0.1 \text{ ng}}{5 \text{ g}} = 0.02 \frac{\text{ng}}{\text{g}} \text{ or } 20 \frac{\text{ng}}{\text{kg}}$$

Sciex OS software is used to prepare a linear standard curve where x is the concentration ratio (analyte/SS) and y is the instrument response ratio (analyte/SS) with 1/x weighting. Surrogates and their internal standard pairs are listed in Table 5, which are used to calculate absolute recoveries of the surrogate standards over the entire extraction method. Surrogates and their native analyte pairs are also listed in Table 5 with their curve fit. The calibration curve has surrogate and internal standard concentrations of 1 ng/mL.

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## 2024.8.4 LC-MS/MS Analysis

All samples were analyzed using a liquid chromatograph (Nexera X2, (Shimadzu, Kyoto, Japan)). The MS/MS data was acquired using scheduled MRM with a Sciex QTRAP® 6500+.


Figure 8-1. Example chromatogram of a spiked single-laboratory validation sample (salmon) with native and labeled PFAS concentrations at 5 µg/kg.

1. PFBA, 2. PFPeA, 3. PFHxA, 4. PFHpA, 5. PFOA, 6. PFNA, 7. PFDA, 8. PFUnA, 9. PFDoA, 10. PFTTrDA, 11. PFTeDA, 12. PFBS, 13. PFPeS, 14. PFHxS, 15. PFHpS, 16. PFOS, 17. PFNS, 18. PFDS, 19. PFUnDS, 20. PFDoS, 21. PFTTrDS, 22. FOSA, 23. DONA, 24. HFPO-DA, 25. 9Cl-PF3ONS, 26. 11Cl-PF3OUdS, 27. 4:2 FTS, 28. 6:2 FTS, 29. 8:2 FTS, 30. 10:2 FTS 31. d5 N-EtFOSAA, 32. M3 PFBA, 33. M3 PFPeA 34. M5PFHxA, 35. M8PFOA, 36. MPFUDa, 37. MPFDoA, 38. M2PFTeDA, 39. M8 FOSA 40. M3 PFBS, 41. M3PFHxS, 42. 13C PFOS, 43. M3 HFPO-DA 44. Labeled 4:2 FTS, 45. Labeled 6:2 FTS, 46. Labeled 8:2 FTS, 47. Labeled 10:2 FTS

## 2024.8.4 LC-MS/MS Analysis

All samples were analyzed using a liquid chromatograph (Nexera X2, (Shimadzu, Kyoto, Japan)). The MS/MS data was acquired using scheduled MRM with a Sciex QTRAP® 6500+.


Figure 8-1. Example chromatogram of a spiked single-laboratory validation sample (salmon) with native and labeled PFAS concentrations at 5 µg/kg.

1. PFBA, 2. PFPeA, 3. PFHxA, 4. PFHpA, 5. PFOA, 6. PFNA, 7. PFDA, 8. PFUnA, 9. PFDoA, 10. PFTTrDA, 11. PFTeDA, 12. PFBS, 13. PFPeS, 14. PFHxS, 15. PFHpS, 16. PFOS, 17. PFNS, 18. PFDS, 19. PFUnDS, 20. PFDoS, 21. PFTTrDS, 22. FOSA, 23. DONA, 24. HFPO-DA, 25. 9Cl-PF3ONS, 26. 11Cl-PF3OUdS, 27. 4:2 FTS, 28. 6:2 FTS, 29. 8:2 FTS, 30. 10:2 FTS 31. d5 N-EtFOSAA, 32. M3 PFBA, 33. M3 PFPeA 34. M5PFHxA, 35. M8PFOA, 36. MPFUDa, 37. MPFDoA, 38. M2PFTeDA, 39. M8 FOSA 40. M3 PFBS, 41. M3PFHxS, 42. 13C PFOS, 43. M3 HFPO-DA 44. Labeled 4:2 FTS, 45. Labeled 6:2 FTS, 46. Labeled 8:2 FTS, 47. Labeled 10:2 FTS

### 2024.12.3 Method Detection Limits

Method detection limits were calculated by performing a minimum of 7 low-level spikes at 0.05 ng/g for all food and feed matrices. The standard deviation of the replicates was multiplied by 3.14 (t-value for seven replicates where  $1-\alpha = 0.99$ ). The MDL is defined as the statistically calculated minimum concentration that can be measured with 99% confidence that the measured concentration is distinguishable from method blank results. This procedure is published in the Code of Federal Regulations.<sup>3</sup> MDLs are recalculated yearly and are dependent on instrumental conditions at the time of the study. The limit of quantification (LOQ) can be calculated by multiplying the standard deviation of the replicates by 10. Both the MDLs and LOQs can be found in Appendix 1D.

| PARAMETER<br>CHỈ TIÊU                                 | METHOD<br>PHƯƠNG PHÁP                                                                                      |
|-------------------------------------------------------|------------------------------------------------------------------------------------------------------------|
| Agricultural chemical residues list of 1110 parameter | LC/MS/MS: Ref. EN 15662:2018<br>GC/MS/MS: Ref. AOAC 2007.01 (21 <sup>st</sup> Ed., 2019 and EN 15662:2018) |

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Selected product: Rice


Main product of the group or subgroup    Other products to which the same MRLs apply

| Code Number | Name | Scientific Name |
|-------------|------|-----------------|
|-------------|------|-----------------|

↓  
 ↓  
 517 chỉ tiêu

↓  
 Maximum residue level Min 0.001 (mg/kg)  
 (2-amino-4-methoxy-6-(trifluormethyl)-  
 1,3,5-triazine (AMTT))

<https://ec.europa.eu/food/plant/pesticides/eu-pesticides-database/start/screen/mrls/latest>

## Quick Method for the Analysis of Highly Polar Pesticides in Food Involving Extraction with Acidified Methanol and LC- or IC-MS/MS Measurement *I. Food of Plant Origin (QuPpe-PO-Method)*

**Version 12.2** (21.12.2023, Document History, see page **104**)

**Check for latest version of this Method under** <https://www.guppe.eu>; older versions: [obsolete versions](#)


**Authors:** M. Anastassiades; A.-K. Schäfer; D. I. Kolberg; E. Eichhorn; H. Dias; A. Benkenstein; S. Zechmann;  
D. Mack; C. Wildgrube; A. Barth; I. Sigalov; S. Goerlich; D. Dörk; G. Cerchia

**EU Reference Laboratory for pesticides requiring Single Residue Methods (EURL-SRM)**

**Contact:** CVUA Stuttgart, Schaflandstr. 3/2, DE-70736 Fellbach, Germany, [www.eurl-pesticides.eu](http://www.eurl-pesticides.eu); [EURL@cvuas.bwl.de](mailto:EURL@cvuas.bwl.de)

**Note:** Changes from V12.1 to V12.2 are highlighted in green

<http://www.eurl-pesticides.eu/>

**Quick Method for the Analysis of  
Numerous Highly Polar Pesticides in Food Involving  
Extraction with Acidified Methanol and LC-MS/MS Measurement  
II. Food of Animal Origin (QuPpe-AO-Method)**

**Version 3.2** (14.05.2019, Document History, see page 24)

Please check for the **latest version of this document** under [www.quppe.eu](http://www.quppe.eu)

Changes from V2 to V3 are highlighted in yellow; from V3 to V3.1 in green and from V3.1 to 3.2 in grey

**Authors:** M. Anastassiades; A.-K. Wachtler; D. I. Kolberg; E. Eichhorn; A. Benkenstein; S. Zechmann;  
D. Mack; A. Barth; C. Wildgrube; I. Sigalov; S. Görlich; D. Dörk; G. Cerchia

**EU Reference Laboratory for pesticides requiring Single Residue Methods (EURL-SRM)**

Address: CVUA Stuttgart, Schaflandstr. 3/2, DE-70736 Fellbach, Germany

Web: [www.eurl-pesticides.eu](http://www.eurl-pesticides.eu),

E-Mail: [EURL@cvuas.bwl.de](mailto:EURL@cvuas.bwl.de)

|                                                                                               |    |
|-----------------------------------------------------------------------------------------------|----|
| 5.6. LC-and IC-MS/MS analysis .....                                                           | 16 |
| 5.6.1. General hints on analytes to avoid pitfalls .....                                      | 22 |
| 5.6.2. Method 1.1 (M1.1): "Gly&Co. AS 11" .....                                               | 33 |
| 5.6.3. Method 1.2 (M1.2): "Gly&Co. AS 11-HC" .....                                            | 35 |
| 5.6.4. Method 1.3 (M1.3): "Gly&Co. Hypercarb" .....                                           | 37 |
| 5.6.5. Method 1.4 (M1.4): "PerChloPhos" .....                                                 | 39 |
| 5.6.6. Method 1.5 (M1.5): "Gly&Co. on Trinity Q1" .....                                       | 41 |
| 5.6.7. Method 1.6 "Gly&Co. on Torus DEA; (M1.6a)" or "Gly&Co. on APPC; (M1.6b)" .....         | 43 |
| 5.6.8. Method 1.7 "PerChloPhos on Torus DEA; (M1.7a)" or "PerChloPhos on APPC; (M1.7b)" ..... | 47 |
| 5.6.9. Method 1.8 (M1.8): "PerChloCyanMalein on Anionic Polar Pesticide Column (APPC)" .....  | 50 |
| 5.6.10. Method 1.9 (M1.9): "Gly&Co. on Raptor Polar X" .....                                  | 51 |
| 5.6.11. Method 1.10 (M1.10): "Gly&Co. on Obelisc N" .....                                     | 53 |
| 5.6.12. Method 1.11 (M1.11): "Gly&Co. on Luna Polar Pesticides" .....                         | 54 |
| 5.6.13. Method 1.12 (M1.12): "PerChloPhos on Luna Polar Pesticides" .....                     | 56 |
| 5.6.14. Method 2 (M2): "Fosetyl and Maleic Hydrazide" .....                                   | 57 |
| 5.6.15. Method 3 (M3): "Amitrole&Co." .....                                                   | 59 |
| 5.6.16. Method 4.1 (M4.1): "Quats&Co. Obelisc R" .....                                        | 61 |
| 5.6.17. Method 4.2 (M4.2): "Quats&Co. BEH Amide" .....                                        | 63 |
| 5.6.18. Method 5 (M5): "Quats&Co. MonoChrom MS" .....                                         | 67 |
| 5.6.19. Method 6 (M6): "Streptomycin and Kasugamycin" .....                                   | 68 |
| 5.6.20. Method 7 (M7): "Morpholine, Diethanolamine and Triethanolamine" .....                 | 69 |
| 5.6.21. Method 8 (M8): "Triazole derivative metabolites (TDMs)" .....                         | 70 |
| 5.6.22. Method 9 (M9): "Difluoroacetic acid and Trifluoroacetic acid" .....                   | 71 |
| 5.6.23. Method 10 (M10): "Triazole derivative metabolites (TDMs) on Torus DEA" .....          | 72 |
| 5.6.24. Method 11 (M11): "Gly&Co. by IC on AS19" .....                                        | 73 |

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EURL-SRM

EU Reference Laboratories for Residues of Pesticides  
Single Residue Methods

✓ = satisfactory chromatography and detection sensitivity achieved,  
NT = Not tested under the conditions shown in the respective sections,  
(✓) = possible but compromised due to matrix effects or lacking separation or limited sensitivity or limitations in the detection of qualifiers compromising identification.  
“-“ analysis was tested and found to be poor under the described conditions,  
\* Using a gradient (98% B -> 60% B in 5 min, hold 2 min)  
\*\* Different LC-conditions required to improve peaks (see M1.7 or M1.9)  
\*\*\* Compromised quantitation of Phosphonic acid due to co-elution of Phosphonic acid and Fosetyl, see also General Hints 5.6.1  
\*\*\*\* Quality of analysis may strongly depend on instrument type and condition

Table 6: Scope of QuPPE-LC-Methods of analytes analyzed in the ESI-neg. mode Part.I

| QuPPE method code                                     | M 1.1<br>(5.6.2) | M 1.2<br>(5.6.3) | M 1.3<br>(5.6.4) | M 1.4<br>(5.6.5) | M 1.5<br>(5.6.6) | M 1.6a/b<br>(5.6.7)      | M 1.7 a/b<br>(5.6.8)     | M 1.8<br>(5.6.9) | M 1.9<br>(5.6.10) | M 1.10<br>(5.6.11) | M 1.11<br>(5.6.12)       | M 1.12<br>(5.6.13)       |
|-------------------------------------------------------|------------------|------------------|------------------|------------------|------------------|--------------------------|--------------------------|------------------|-------------------|--------------------|--------------------------|--------------------------|
| Separation principle                                  | Anion Ex.        | Anion Ex.        | Carbon           | Carbon           | HILIC            | HILIC                    | HILIC                    | HILIC            | HILIC             | HILIC              | HILIC                    | HILIC                    |
| Column type                                           | AS 11            | AS 11-HC         | Hypercarb        | Hypercarb        | Trinity Q1       | a: Torus DEA/<br>b: APPC | a: Torus DEA/<br>b: APPC | APPC             | Raptor<br>PolarX  | ObeliscN           | Luna Polar<br>Pesticides | Luna Polar<br>Pesticides |
| ANALYTES COVERED BY LC-MS/MS IN THE ESI-NEGATIVE MODE |                  |                  |                  |                  |                  |                          |                          |                  |                   |                    |                          |                          |
| Ethephon                                              | ✓                | ✓                | ✓                | NT               | ✓                | ✓                        | (✓)                      | NT               | ✓                 | ✓                  | ✓                        | -                        |
| HEPA                                                  | ✓                | ✓                | ✓                | NT               | ✓                | ✓                        | (✓)                      | NT               | ✓                 | ✓                  | ✓                        | -                        |
| Glufosinate                                           | ✓                | ✓                | ✓                | NT               | ✓                | ✓                        | (✓)                      | NT               | ✓                 | ✓                  | ✓                        | -                        |
| N-Acetyl-Glufosinate                                  | ✓                | ✓                | ✓                | NT               | ✓                | ✓                        | (✓)                      | NT               | ✓                 | ✓                  | ✓                        | -                        |
| MPPA                                                  | ✓                | ✓                | ✓                | NT               | ✓                | ✓                        | (✓)                      | NT               | ✓                 | ✓                  | ✓                        | -                        |
| Glyphosate                                            | ✓                | ✓                | ✓                | NT               | ✓                | ✓                        | (✓)                      | NT               | ✓                 | ✓                  | ✓                        | -                        |
| AMPA                                                  | ✓                | ✓                | ✓                | NT               | ✓                | ✓                        | (✓)                      | NT               | ✓                 | (✓)                | (✓)                      | -                        |
| Phosphonic acid                                       | (✓)              | (✓)              | ✓                | ✓                | ✓                | ✓                        | ✓                        | NT               | ✓                 | (✓)***             | (✓)                      | ✓                        |
| N-Acetyl-AMPA                                         | NT               | ✓                | ✓                | NT               | ✓                | ✓                        | (✓)                      | NT               | NT                | NT                 | NT                       | NT                       |
| Fosetyl-Al                                            | -                | ✓                | ✓                | NT               | ✓                | ✓                        | (✓)                      | NT               | ✓                 | (✓)***             | ✓                        | -                        |
| Maleic Hydrazide                                      | -                | -                | ✓                | NT               | -                | -                        | -                        | (✓)              | (✓)               | (✓)                | -                        | -                        |
| Perchlorate                                           | NT               | -                | ✓                | ✓                | ✓                | (✓)**                    | ✓                        | ✓                | ✓                 | ✓                  | -                        | ✓                        |
| Chlorate                                              | NT               | -                | ✓                | ✓                | ✓                | (✓)**                    | ✓                        | ✓                | (✓)               | ✓                  | -                        | (✓)                      |
| Bialaphos                                             | NT               | NT               | ✓                | NT               | ✓                | NT                       | NT                       | NT               | ✓                 | NT                 | NT                       | NT                       |
| Cyanuric acid                                         | NT               | NT               | ✓                | NT               | -                | -                        | -                        | ✓                | (✓)               | (✓)                | -                        | -                        |
| Bromide                                               | NT               | NT               | -                | ✓                | ✓                | (✓)**                    | ✓                        | NT               | ✓                 | ✓                  | (✓)                      | ✓                        |
| Bromate                                               | NT               | NT               | (✓)              | ✓                | NT               | NT                       | ✓                        | NT               | NT                | NT                 | NT                       | NT                       |
| N-Acetyl-Glyphosate                                   | NT               | NT               | ✓                | NT               | (✓)**            | ✓                        | (✓)                      | NT               | -                 | ✓                  | ✓                        | -                        |
| Difluoroacetic acid                                   | NT               | NT               | NT               | NT               | NT               | ✓                        | ✓                        | NT               | ✓                 | NT                 | NT                       | NT                       |
| Trifluoroacetic acid                                  | NT               | NT               | NT               | NT               | NT               | ✓                        | ✓                        | NT               | ✓                 | ✓                  | -                        | (✓)                      |
| Thiocyanate                                           | NT               | NT               | (✓)**            | ✓                | NT               | NT                       | ✓                        | NT               | ✓                 | NT                 | -                        | (✓)                      |
| Desmethyl-Dimethoate                                  | NT               | NT               | ✓                | -                | NT               | NT                       | NT                       | NT               | NT                | NT                 | NT                       | NT                       |
| Diethyl thiophosphate                                 | NT               | NT               | NT               | NT               | NT               | NT                       | ✓                        | NT               | NT                | NT                 | NT                       | NT                       |
| Diethyl phosphate                                     | NT               | NT               | NT               | NT               | NT               | NT                       | ✓                        | NT               | NT                | NT                 | NT                       | NT                       |
| O,O-dimethyl dithiophosphate                          | NT               | NT               | NT               | NT               | NT               | NT                       | ✓                        | NT               | NT                | NT                 | NT                       | NT                       |
| O,O-dimethyl thiophosphate                            | NT               | NT               | NT               | NT               | NT               | NT                       | ✓                        | NT               | NT                | NT                 | NT                       | NT                       |
| Diethyl dithiophosphate                               | NT               | NT               | NT               | NT               | NT               | NT                       | ✓                        | NT               | NT                | NT                 | NT                       | NT                       |
| Dimethyl phosphate                                    | NT               | NT               | NT               | NT               | NT               | NT                       | ✓                        | NT               | NT                | NT                 | NT                       | NT                       |

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

## EURL-SRM EU Reference Laboratories for Residues of Pesticides Single Residue Methods

**Table 8 : Methods mainly used by CVUA Stuttgart**

| Method                                                      | Special remarks on Substances                                                                                                                                                                                                                                                                            | LC-MS/MS                                       | Comments                                                                                                                                                                                            |
|-------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>M 1.3: Glyphosate &amp; Co. Hypercarb (see 5.6.4)</b>    | Glyphosate<br>AMPA<br>N-Acetyl-AMPA<br>N-Acetyl-Glyphosate<br>Ethephon<br>HEPA<br>Glufosinate<br>N-Acetyl-Glufosinate<br>MPPA<br>Fosetyl-Al<br>Phosphonic acid (screening)<br>Maleic Hydrazide<br>Perchlorate (screening))<br>Chlorate (screening)<br>Cyanuric acid<br>Bialaphos<br>Desmethyl-Dimethoate | 1290 Agilent Infinity II and Sciex QTRAP 6500+ | Evaluation via solvent calibration and IL-ISs except for Bialaphos and N-Acetyl-AMPA (IL-IS not yet available)<br>M 1.5 and M 1.6 are currently being tested for their suitability to replace M 1.3 |
| <b>M 1.4: Method 1.4 (M1.4): "Per-ChloPhos" (see 5.6.5)</b> | Perchlorate (quantitative)<br>Chlorate (quantitative)<br>Phosphonic acid (quantitative)<br>Bromide (Screening, quantitative)<br>Bromate (quantitative)<br>Thiocyanate                                                                                                                                    | 1290 Agilent Infinity II and Sciex QTRAP 6500+ | Mostly employed directly (option: screening by M 1.3, if positive -> M 1.4)<br>Dilution 5-fold<br>Evaluation via solvent calibration and IL-ISs                                                     |

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**Figure 9: Typical chromatograms of Ethephon, HEPA, Glyphosat, AMPA, Glufosinate, MPPA, N-Acetyl-AMPA, N-Acetyl-Glufosinate, Fosetyl-Al and Phosphonic acid at 0.1 mg/L in methanol with 1% formic acid.**

**Oat**  
0.034 ppm

**Figs**  
0.12 ppm

**Pineapple**  
0.009 ppm

## Biggest Problem with LC- and GC- MS(/MS)


# Matrix Effects



plus costs to purchase and maintain,  
and facility requirements, and downtime,  
and need for more expertise due to  
greater complexity


It is still a pain!

## How to Calculate (Estimate) Matrix Effects


[https://www.researchgate.net/figure/Matrix-effect-for-the-six-pesticides-studied-in-the-bell-pepper-matrix\\_fig2\\_335403723](https://www.researchgate.net/figure/Matrix-effect-for-the-six-pesticides-studied-in-the-bell-pepper-matrix_fig2_335403723)


## Effect of Analyte Protectants


## Analyte Protectants

Strongly interact with active sites in GC system (inlet, column and ion source) to decrease degradation and adsorption of co-injected analytes.


Sharper peaks, less tailing, more ruggedness, lower LOD


**ethylglycerol**  
1 mg/mL


**gulonolactone**  
0.1 mg/mL

**sorbitol**  
0.1 mg/mL

**shikimic acid**  
0.05 mg/mL

**Fig 2.**

Optimization of analyte protectant (AP) for GC-MS/MS signal enhancement of dimethipin after priming with pepper leaf matrix (a) no priming (b) 10% AP after priming (c) 20% AP after priming (d) 30% AP after priming.

## SCION 8900 TQMS APPLICATION NOTE

### Multiplexed pesticide analysis by GC-TQMS ANALYSIS OF 450 PESTICIDES IN BLUEBERRIES


Figure 3. Increasing MRM does not impede sensitivity

Despite the increase in number of components measured, (Figure 3), the SCION 8900 TQMS maintains good accuracy and linearity at scan times of 4ms. The ion ratio stability is maintained when switched from 40ms to 7ms scan time. This is maintained even at the lowest concentration of calibration standards (down to 1ppb). Reproducibility is maintained over multiple runs.

minute and held for 14.3 minutes. The Backflush option was activated with a pressure of 5PSI and a start time of 14 minutes.

### RESULTS

The data showed excellent linearity for standards and QC samples (Figure 5) with perfect ion scores observed throughout the batch. Retention time reproducibility was extremely consistent and good recovery and reproducibility. This reproducibility was consistent through all of the detected analytes. For example, chlorpheninfos, an organophosphorous compound that is widely used as an insecticide, was detected at a concentration of 8ppb. The retention time of the compound was 8.45 minutes, which is within 0.02 minutes of the expected RT. The calibration curve had a  $R^2$  of >0.998, with a residual of less than  $\pm$  20%.

| Name                | RT   | RT   | RT   | RT   | RT   | RT   | RT   | RT   | RT   |
|---------------------|------|------|------|------|------|------|------|------|------|
| 1. Acetone          | 1.1  | 1.1  | 1.1  | 1.1  | 1.1  | 1.1  | 1.1  | 1.1  | 1.1  |
| 2. Ethanol          | 1.2  | 1.2  | 1.2  | 1.2  | 1.2  | 1.2  | 1.2  | 1.2  | 1.2  |
| 3. Chloropheninfos  | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 |
| 4. Chloropheninfos  | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 |
| 5. Chloropheninfos  | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 |
| 6. Chloropheninfos  | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 |
| 7. Chloropheninfos  | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 |
| 8. Chloropheninfos  | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 |
| 9. Chloropheninfos  | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 |
| 10. Chloropheninfos | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 | 8.45 |

Direct aqueous analysis of pesticides and PPCPs in drinking and bottled water at parts per trillion levels

## Samples

- Pesticides and PPCP in water
- Compounds: positive mode 370, negative mode 61 (46 internal stds)
- The calibration standard concentrations were 0.1, 0.2, 0.5, 1, 2, 5, 10, 20, 50, 100, 200 and 500 ng/L. A calibration blank was also prepared.
- MilliQ water, tap water (drinking water), Evian, FIJI

## Methods

LC system: ExionLC™ System

- Column: Luna Omega C18 1.6  $\mu\text{m}$ , 100 $\times$ 2.1 mm
- Mobile phase:
  - A : Water + 0.1% formic acid + 5mM NH<sub>4</sub> formate
  - B : MeOH + 0.1 % formic acid + 5mM NH<sub>4</sub> formate
- Temperature: 40 °C
- Injection volume: 500  $\mu\text{L}$
- MS: fast pos/neg switching with Scheduled MRM™ Algorithm

| Total Time (min) | Flow Rate ( $\mu\text{L}/\text{min}$ ) | A (%) | B(%) |
|------------------|----------------------------------------|-------|------|
| 0.00             | 500                                    | 99    | 1    |
| 0.50             | 500                                    | 99    | 1    |
| 16.00            | 500                                    | 1     | 99   |
| 21.00            | 500                                    | 1     | 99   |
| 21.10            | 500                                    | 99    | 1    |
| 25.00            | 500                                    | 99    | 1    |


## SCIEX QTRAP®® 6500+

| Product                                   | Glufosinate sum |     |             | Fosetyl sum |      |             | Glyphosate |     |             | Chlorate |     |             | Ethephon |     |             |
|-------------------------------------------|-----------------|-----|-------------|-------------|------|-------------|------------|-----|-------------|----------|-----|-------------|----------|-----|-------------|
|                                           | LOD             | MRL | %RSD at MRL | LOD         | MRL  | %RSD at MRL | LOD        | MRL | %RSD at MRL | LOD      | MRL | %RSD at MRL | LOD      | MRL | %RSD at MRL |
| <i>Fruit and Vegetables</i>               | 16              | 30  | 11%         | 25          | 2000 | 13%         | 5          | 100 | 15%         | 8        | 10  | 15%         | 18       | 50  | 11%         |
| <i>Seeds</i>                              | 12              | 30  | 12%         | 90          | 2000 | 15%         | 8          | 100 | 15%         | 3        | 10  | 10%         | 6        | 50  | 14%         |
| <i>Vegetable oil, Fat and Fatty Acids</i> | 15              | 30  | 19%         | 40          | 2000 | 12%         | 7          | 100 | 22%         | 2        | 10  | 6%          | 3        | 50  | 7%          |
| <i>Grain</i>                              | 18              | 30  | 12%         | 71          | 2000 | 14%         | 8          | 100 | 7%          | 7        | 10  | 14%         | 9        | 50  | 6%          |
| <i>Herbs and spices</i>                   | 25              | 100 | 8%          | 87          | 2000 | 13%         | 23         | 100 | 6%          | 8        | 10  | 15%         | 8        | 100 | 16%         |
| <i>Meat and Seafood</i>                   | 19              | 30  | 15%         | 23          | 100  | 12%         | 9          | 50  | 23%         | 4        | 10  | 8%          | 4        | 50  | 10%         |
| <i>Animal Oil, Fat and Fatty Acids</i>    | 14              | 30  | 20%         | 51          | 100  | 11%         | 9          | 50  | 25%         | 10       | 10  | 16%         | 7        | 50  | 12%         |
| <i>Eggs and Eggs products</i>             | 18              | 30  | 12%         | 33          | 100  | 11%         | 4          | 50  | 13%         | 12       | 10  | 9%          | 6        | 50  | 17%         |
| <i>Milk and Milk products</i>             | 17              | 30  | 9%          | 20          | 100  | 6%          | 8          | 50  | 22%         | 5        | 10  | 12%         | 5        | 50  | 13%         |
| <i>Fatty acids</i>                        | 21              | 100 | 14%         | 70          | 1000 | 14%         | 3          | 100 | 18%         | 4        | 10  | 9%          | 3        | 100 | 10%         |


- For Fast Polarity Switching **IonDrive High Energy Dynode detector – 6500+**

## IonDrive HED detector

- Positive Ions strike the HED first, Negative ions strike the magnum directly
- Magnum to be kept at one polarity


Positive ions first converted to an electron, then electrons are detected by electron multiplier



Negative ions directly detected by electron multiplier

## CEM detector

- Both positive and negative ions strike to detector
- Detector voltage must be switched from positive to negative through ground for each cycle

## TAKE ADVANTAGE BY SCANNING MORE COMPOUNDS

Gain of 90 msec can be used to add more compounds to each scan, or dwell longer on existing MRMs to improve ion statistics

| PARAMETER<br>CHỈ TIÊU                                 | METHOD<br>PHƯƠNG PHÁP                                                                                      |
|-------------------------------------------------------|------------------------------------------------------------------------------------------------------------|
| Agricultural chemical residues list of 1110 parameter | LC/MS/MS: Ref. EN 15662:2018<br>GC/MS/MS: Ref. AOAC 2007.01 (21 <sup>st</sup> Ed., 2019 and EN 15662:2018) |

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
XIC of positive and negative ionization

SCIEX 7500 System

ALL QUANTIFIER IONS IN MILLI-Q WATER

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For Research Use Only. Not for Use In Diagnostic Procedures.

## Food and Environmental

### Direct analysis of polar pesticides in water: a routine accredited method

*Using the SCIEX QTRAP® 6500+ LC-MS/MS System*

Y. Fillâtre<sup>1</sup>, J. Prades<sup>2</sup>, L. Ley<sup>2</sup>, D. McMillan<sup>1</sup>

<sup>1</sup>SCIEX, France, <sup>2</sup>LDAR24, France

"The Glyphosate Paradox" is that, although glyphosate is the most widely used agrochemical in the world, it is also one of the least often determined by analytical methods.<sup>1</sup> Monitoring a highly polar, small organic pesticide such as glyphosate in food and water from diverse sources presents a significant challenge. Polar pesticides are not amenable to standard extraction procedures, are frequently poorly ionized, and demonstrate poor chromatographic separation. These pesticides, therefore, have historically required complex, single-residue methods to make them amenable to analysis—usually involving time-consuming derivatization steps and considerable clean-up procedures.

In 2018, Nefelab, in collaboration with SCIEX, developed a

and other polar pesticides in water samples. In France, to become accredited by the COFRAC (Comité Français d'Accréditation), the laboratory must demonstrate their analytical method is fit for a routine use by achieving required

### Key features


- Accredited method in drinking water with  
✓ LOQ = 20 ng/L for fosetyl-Al  
✓ LOQ = 30 ng/L for AMPA, glufosinate, glyphosate
- Flexibility and ease of use
- Long term stability and robustness

**Methods Sample preparation:** Chlorinated water samples were stabilized with sodium thiosulfate to neutralize the presence of chlorine and then analyzed with the LC-MS/MS system.

**Surface water and underground water were directly injected.**

**Chromatography:** Separation was achieved using an ExionLC™ System fitted with a 500 µL loop and a CTO-20A column oven.

A LC column of 150 x 4 mm with a guard column of the same material and a 0.5 µm filter was used for separation. This column allows the use of MS amenable mobile phases at around pH 9.

**Figure 3. Calibration curves obtained for the third accreditation sequence.** (Top) Calibration curves of quantification and qualification MRMs of AMPA, fosetyl-Al, glufosinate and glyphosate corrected with internal standards for S3. (Bottom) Extracted ion chromatogram (XIC) of MRM 1 and MRM 2 for the blank, 20 ng/L and 30 ng/L samples for the four compounds. The solid blue line shows the mean ion ratio calculated from the standards and the dotted lines the tolerance of 20%.

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^^^ SCIEX OS  
^^^ Software A SIMPLE SWITCH  
^^^ FOR A BIG DIFFERENCE

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## Drug Discovery and Development

# Quantitative LC-MS Solution for Targeted Analysis of Cell Culture Media

*Featuring the SCIEX QTRAP® 6500+ LC-MS/MS System*

Zuzana Demianova<sup>1</sup> and Lei Xiong<sup>2</sup>

<sup>1</sup>SCIEX, Brea, CA; <sup>2</sup>SCIEX, Redwood City, CA

Measurement of polar and nonpolar as well as positive and negative polarity components within a single method

| Component group | Number of components |
|-----------------|----------------------|
| Amino acids     | 39                   |
| Vitamins        | 15                   |
| Carbohydrates   | 4                    |
| Fatty acids     | 5                    |
| Nucleic acids   | 17                   |
| Others          | 32                   |

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## Drug Discovery and Development

# Quantitative LC-MS Solution for Targeted Analysis of Cell Culture Media

*Featuring the SCIEX QTRAP® 6500+ LC-MS/MS System*

Zuzana Demianova<sup>1</sup> and Lei Xiong<sup>2</sup>

<sup>1</sup>SCIEX, Brea, CA; <sup>2</sup>SCIEX, Redwood City, CA

Measurement of polar and nonpolar as well as positive and negative polarity components within a single method

| Component group | Number of components |
|-----------------|----------------------|
| Amino acids     | 39                   |
| Vitamins        | 15                   |
| Carbohydrates   | 4                    |
| Fatty acids     | 5                    |
| Nucleic acids   | 17                   |
| Others          | 32                   |

**Mass Spec:** MRM parameters for 110 components

- + a Triple Quad™ 6500+ system.
- + The scheduled MRM™ Algorithm Pro(sMRM) was used

- + Two MRM transitions were monitored for each analyte,
- + This method contains 178 positive MRMs and 54 negative MRMs.

**Chromatography:** an ExionLC™ System

- + Phenomenex Kinetex® F5 column (150 mm x 2.1 mm ID, particle size 2.6 µm).
- + Total method time was 20 min
- + flowrate of 200 µl/min.

Mobile phase A was 0.1% formic acid in water, and mobile phase B was acetonitrile with 0.1% formic acid.

- + Column oven temperature was 40°C.
- + Injection volume was 5 µL.

**Sample preparation:** A final master mix was prepared containing all standards (6.67-20 µg/ml) depending on analyte signal).

Cell culture medium (CD CHO medium, Gibco) was diluted 5-fold with 0.1% formic acid in 50% acetonitrile, and centrifuged. The supernatant was further diluted 60-fold with 0.1% formic acid.


Figure 5. Representative MIDAS™ workflow for L-Leucine and L-Isoleucine with compound match from the Sciex Accurate Metabolite Spectra

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


Figure 3. Top panel shows various XICs of low abundant components from Gibco's cell culture media (CD CHO medium). Bottom panel shows amino acids from Gibco's cell culture media (CD CHO medium).

High sensitivity  
Selectivity  
Speed  
Robustness;  
Enables unambiguous  
identification and  
quantification of a large  
number of analytes in a  
single analysis

## SCIEX CCMA Method Contains > 110 Media Analytes

## Food and Environmental

**Table 1.** Daily required values (DV) of different B vitamins for a human adult as obtained from the FDA<sup>1</sup>

| Vitamin          |     | DV (mg) |
|------------------|-----|---------|
| Thiamine         | B1  | 1.5     |
| Riboflavin       | B2  | 1.7     |
| Niacin           | B3  | 20      |
| Pantothenic acid | B5  | 10      |
| Pyridoxal        | B6  | 2       |
| Biotin           | B7  | 0.3     |
| Folic acid       | B9  | 0.4     |
| Cyanocobalamin   | B12 | 0.006   |

## Complex in Infant Formula

Concord, Ontario (Canada)


Vitamin B is a group of water-soluble vitamins

Food samples: matrices complex

Sensitive methods typically require highly selective sample clean up procedures.

Vitamin B is a complex mixture of highly polar compounds (Figure 1) whose  $pK_a$  range from 0.5 to 10.2, making their analysis challenging.

**Figure 1.** Chemical structures of B vitamins

## Sample Preparation

Sample (1 g) was mixed with 50% acetonitrile in acidified water (containing an antioxidant) + IS.

Shaken vigorously for 1 minute, roller mixed for 10 minutes (protected from light).

After centrifugation, filtered and the filtrate diluted 20 times with water containing an ion paring reagent. The sample preparation was kept as simple as possible to reduce possible vitamin breakdown, with SPE no longer needed for the late eluting B7,

B9, and B12 vitamins.

During the development work the effects of light, temperature, and acidity on standard stability were tested and it was found that the use of amber glass with a lower pH with the presence of an antioxidant helped stabilize the extracts.