



Detection of fluoride and chloride in germanium ores

Introduction:

Germanium has various special properties and has extensive and important applications in fields such as semiconductors, aerospace measurement and control, nuclear physics detection, fiber optic communication, infrared optics, solar cells, chemical catalysts, biomedicine, etc. It is an important strategic resource. In the electronic industry, in alloy pretreatment, and in the optical industry, it can also serve as a catalyst.

Germanium products will further increase technological innovation and transformation efforts, optimize process flow, and improve resource recovery and utilization rates. Due to its extensive demand and industrial upgrading and transformation, there is currently a demand for quality testing of germanium ore. The content of fluorine and chlorine in germanium ore is an important indicator, as fluorine and chlorine not only pollute the environment, but also corrode equipment and pipelines during the refining process, affecting product quality. Therefore, the research on the detection and methods of fluoride and chloride content in germanium ore has attracted widespread attention.

Detection items (Table 1):

| | | |
|--------------|----------------|-----------------|
| Anion | F ⁻ | Cl ⁻ |
|--------------|----------------|-----------------|

Keywords: On-line ion chromatography, Fluoride, Chloride

Instruments and equipment

- **Ion chromatograph:** SH-CIC-3200
- **Ultra pure water machine:** ECO-S15


Qingdao Shenghan Chromatograph Technology Co., Ltd

Requirements

Reagents

Unless otherwise specified, all reagents used are superior grade. F⁻ Cl⁻ anions standard solution (1000 mg/L)

Deionized Water



More information, Please visit our website:
<http://www.sheng-han.net/>
 Serial number:051

When preparing standard samples manually or diluting real samples, please use ASTM filtration and deionization requirements that meet the specifications listed in the table 2.

Table 2: Deionized water specification.

| Specification | |
|------------------------|-------------|
| Ions Resistivity | ≥18.25MΩ·cm |
| Organics-TOC | <10ppb |
| Iron/Transition Metals | <1ppb |
| Pyrogens | <0.03Eu/mL |
| Particulates (>0.2μm) | <1unit/mL |
| Colloids-Silica | <10ppb |
| Bacteria | <1cfu/mL |

Chromatography conditions (Anions):

Table 3: Anions analysis conditions

| | |
|-------------------------------|-------------|
| Instrument | SH-CIC-3200 |
| Eluent | 15 mM KOH |
| Flow rate | 1.0 mL/min |
| Injection volume | 25 μL |
| Analytical Column | SH-AC-23 |
| Column oven temperature | 35°C |
| Conductivity cell temperature | 35°C |
| Suppressor current | 45 mA |

Sample preparation

First, grind the sample to a particle size not greater than 0.080 mm, and dry the sieved sample in a 105 °C ± 2 °C oven for 1 hour. After removal, cool it to room temperature in a dryer for later use.

Table 4: Sample preparation

| No. | Weight (g) | Volume (mL) | Diluting solvent |
|-----|------------|-------------|------------------|
| 1# | 0.0998 | 100 | NaOH solution |
| 2# | 0.0997 | 100 | NaOH solution |

Standard chromatogram

Standard chromatogram, As shown in below:




Figure 1. Chromatogram of standard sample.

Blank chromatogram




Figure 2. Chromatogram of blank

Sample chromatogram




Figure 3. Chromatogram of sample 1#




Figure 4. Chromatogram of sample 2#

Results and calculations

Table 4: Sample test result

| Sample | F ⁻ (mg/L) | Cl ⁻ (mg/L) |
|--------|-----------------------|------------------------|
| 1# | 474.7 | 64.23 |
| 2# | 1780 | 1972 |

Remarks: ① The measured value has been deducted from the blank value; ② There may be differences in testing results between different methods and laboratories.

Feasibility analysis and conclusion

The above experiments prove that the detection method has good resolution and is suitable for the determination of the content of the components to be measured in the sample.