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Model: **IC-8100ST**

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
High-Speed Ion Chromatography Features of IC-8100 Series

High-throughput analysis within 5 minutes

In combination with a high-speed separation column, measurement of anions and cations can be completed within 5 minutes.


- The analysis time can be greatly shortened.
- The measurement time of anions can be shortened from 16 minutes to 5 minutes (compared to our conventional product).
- The measurement time of cations can be shortened from 20 minutes to 5 minutes (compared to our conventional product).

Example of Standard Anions Measurement

Flow rate: Support 0.8 mL/min (Downer) 1.5 mL/min
Pressure: Support 0.5 MPa (Downer) 10.0 MPa

Example of Standard Cations Measurement


Flow rate: Support 0.7 mL/min (Downer) 1.2 mL/min
Pressure: Support 0.9 MPa (Downer) 17.8 MPa

Eluent and Flow Gradient function

Eluent step gradient and flow step gradient are possible. It is convenient for simultaneous measurement of wide range retention capacity ions. And the cleaning of column can be efficiently performed.

- Analysis time is greatly shortened by eluent step gradient.
- Column cleaning can be efficiently performed by flow step gradient.
- The degasser has two flow lines.

Analysis under eluent gradient condition


Column: TSKgel SuperIC-Anion H5 (4.6 mm ID x 10 cm)
TSKgel SuperIC-Cation H5 (4.6 mm ID x 10 cm)
Eluent A: 0.1 M NaOH, 0.5 mmol/L NaCO₃
Eluent B: 0.1 M NaOH, 0.5 mmol/L NaCO₃
Flow rate: 1.0 mL/min (Downer)
Suppressor charger: 0.5 min, 36.5 min
Suppressor gel: TSKgel suppress IC-A
Detection: Conductivity
Injection volume: 20 µL


Automatic gel replacing suppressor system

Stable and continuous high sensitivity analysis can be achieved by gel suppressor system which consists of suppressor gel and rotary valve.

- Suppressor gel is automatically replaced according to the injection timing of autosampler.
- Regeneration of suppressor is not needed.
- Contaminants from sample do not accumulate in suppressor.
- There is no decrease in sensitivity due to deterioration of the suppressor.
- Low dead volume reduces sample diffusion.
- Suppressor gel low sensor is equipped.
- Two types of suppressor gel bottles are available. (30 mL or 60 mL)


Flow of Suppressor Gel Replacement

Suppressor Gel Bottle (60 mL)

Features of Gel Suppressor System

The gel suppressor system automatically conducts the exchange of suppressor gel according to the injection timing of autosampler. Since the suppressor gel is changed after each measurement, it is not affected by sample contaminants that adsorbed on the suppressor gel. Therefore, it is possible to perform analysis at a stable baseline every time.

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
Streamlining routine analysis and variety of functions IC-8100 series dedicated program pursuing "ease of use"

IC-8100-WS workstation is a program with system control, data analysis, and data management functions. Two systems can be controlled, analyzed, and managed with 1 PC.

Easy operation with 3 applications and various tools

- User manager application / user authority setting and registration management
- Acquisition application / system control, data acquisition
- Analysis application / Data analysis, data management
- Log management tool / log confirmation, audit trail tracking
- Validation tool / System evaluation such as noise, detection limit, trueness, precision etc.

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