

LC 6000

Specification sheet

SCION 6310/6320 COLUMN OVEN

Dimensions and Weights

Dimensions: 490.5W x 360D x 174H (mm)

Weight : 13 kg

Power Supply: AC100 to 240 V (50/60 Hz)
(230 VA when optional switching valve is loaded)

Temperature

Control System: 6310: Air circulating system
6320: Peltier heating/cooling system

Setting Range: 1°C to 85°C (0.1°C step)

Accuracy: ±0.5°C

Stability: ±0.1°C

Safety Feature: Thermal protector (105°C) –
Overheat prevention mechanism

Operating temperature range: 4°C to 35°C



Column Storage Capacity

300 mm x 3 (maximum) with guard column

Liquid Leak Sensor

System will alarm if a leak is detected

SCION Instruments HQ

Suite 34, Livingston Business Centre
Kirkton Road South, Livingston,
EH54 7FA, Scotland, UK.
Tel: +44 1506 300 200
sales-eu@scioninstruments.com
www.scioninstruments.com

SCION Instruments NL BV

Amundsenweg 22-24,
4462 GP Goes,
The Netherlands
Tel: +31 (0) 113-287 600
sales-eu@scioninstruments.com
www.scioninstruments.com