

# LC 6000

Specification sheet

## SCION 6510 ORGANISER (SOLVENT CABINET)

### Dimensions & Weights

Dimensions: 340W x 420D x 200H (mm)

Weight : 9 kg

### Output Power

One pump

One autosampler

Two detectors (UV, UV-Vis, Diode Array or RI)


One interface control board

Power Supply/ Consumption: AC 100 V to 240 V (50/60 Hz), 520 VA

### Bottle Capacity

Extra-large solvent cabinet

6 x 1 L bottles and 3 x 500 mL bottle



#### SCION Instruments HQ

Suite 34, Livingston Business Centre  
Kirkton Road South, Livingston,  
EH54 7FA, Scotland, UK.  
Tel: +44 1506 300 200  
sales-eu@scioninstruments.com  
www.scioninstruments.com

#### SCION Instruments NL BV

Amundsenweg 22-24,  
4462 GP Goes,  
The Netherlands  
Tel: +31 (0) 113-287 600  
sales-eu@scioninstruments.com  
www.scioninstruments.com