




# Aquisys 2008

21 CFR Part 11  
compliant Software






Main screen



Measuring menu



Access rights menu

Aquisys 2008 is the ideal supplement for the SCHMIDT+HAENSCH measuring instruments. This software is completely 21 CFR Part 11 compliant and enables you to control your Polarimeters, Refractometers and Density Meters with your PC.

- The integrated, though independent program VALIDATE checks the integrity of Aquisys 2008 before run
- Remote operation of SCHMIDT+HAENSCH Polarimeters, Refractometers and Density Meters with your PC
- Single and repeated measurements
- Audit trail function (history logging, one or all users divided by event type)
- User management with individually configurable access rights for every user
- Logging of electronic recordings and modifications including user information as well as date and time stamp
- Report generation
- Electronic signature which is unmistakably attributed to a person
- The system is able to generate complete paper prints of the electronic records
- Encrypted electronic data storage and archiving
- Backup function (complete settings, audit trail and reports)

#### Typical applications of the model

- Mostly used in pharmaceutical or large laboratories
- In laboratories with plenty of instruments and users
- For quality control of production

