

## WD7015

---



15 DIN BASKETS DOUBLE DOOR WASHER DISINFECTOR  
DOUBLE DOOR WASHER DISINFECTOR FOR CSSD  
DOWNLOAD AREA

---

- ✓ washing up to 95°C - antibacterial thermal disinfection
- ✓ capacity of 18 DIN baskets
- ✓ programs: 20 default + 20 custom
- ✓ multi-language TFT colour touch screen on both sides
- ✓ functioning protected by different levels password
- ✓ integrated built-in panel printer (optional)
- ✓ USB port for downloading washing cycles log
- ✓ n. 2 motors forced hot air drying system
- ✓ Ethernet port for LAN connection (optional)
- ✓ n. 2 sliding full glass doors
- ✓ doors interlocked for complete separation dirty/clean side
- ✓ n. 2 spray arms, one on the top and on the bottom of the chamber
- ✓ high pressure hot air drying system with double absolute 99.999% HEPA filters
- ✓ high efficiency steam condenser (optional)
- ✓ peristaltic pump for liquid detergent and level sensor
- ✓ peristaltic pump for liquid neutralizer and level sensor
- ✓ possibility of installing up to 3 additional peristaltic pumps
- ✓ detergent cabinet with removable tray for storing up to 4 cans of 5lt and one of 1lt
- ✓ washing chamber made of AISI 316L stainless steel
- ✓ external panels made of AISI 304 stainless steel
- ✓ external dimensions LxDxH: 900 x 1000 x 1942 mm
- ✓ internal working dimensions LxDxH: 626 x 812 x 685 mm
- ✓ cold water, hot water and demineralized water connections
- ✓ n. 2 independent washing pumps

- ✓ floor drain connection
- ✓ drain pump (optional)
- ✓ electrical connection: 3/N/PE 380-400 V - 50 Hz - 20 kW max
- ✓ 60 Hz version available
- ✓ without trolleys

Smeg S.p.A. Via Leonardo Da Vinci, 4 - 42016 Guastalla (RE) Italy - P. IVA 01555030350  
Phone +39 0522 8211 - fax +39 0522 821453 - [service.instruments@smeg.it](mailto:service.instruments@smeg.it)

