

## OVERVIEW

The ZA4000 Series Polarized Zeeman Atomic Absorption Spectrophotometer is now available as a successor to the ZA3000 Series Atomic Absorption Spectrophotometer that has earned a high evaluation from the customers.

The ZA4000 Series is available in three different models --- tandem version (ZA4000), flame version (ZA4300), furnace version (ZA4700) and rapid sequential flame (ZA4800). All three versions utilize the polarized Zeeman dual detection system, and feature high sensitivity and precision over a wide range of analyses.

The ZA4000 series Polarized Zeeman Atomic Absorption Spectrophotometer preserve the basic performance of the high precision and sensitive analysis of ZA3000 series, the series incorporates the flame sequential technology as the new technology.

The ZA4000 Series is the new ultimate analysis equipment which considered the improvement of analytic reliability and function, and the environmental performance.



The ZA4000 Series Atomic Absorption Spectrophotometers